

Getinge Clean

Manual Detergent High Foam

Manual Detergent High Foam is a highly effective liquid formula which contains a powerful blend of anionic and nonionic detergents. It quickly penetrates organic soils, emulsifies blood and protein, and is free rinsing for film-free results. Manual Detergent High Foam noncorrosive and safe for soft metals.



Application

For manual prewashing of surgical instruments or equipment. Not for machine use.

Physical and Chemical Properties

- Clear pink, viscous liquid
- Mild, floral fragrance
- pH: 5.0 to 9.5
- Store below 120°F (49°C)

NOTICE: For institutional and industrial use only. Not for home use. Refer to SDS for more information regarding safety, handling and storage and waste disposal methods.

Expiration/Use-By Date

23 Months from date of manufacture.

Key Features

- Neutral pH, mild to skin
- Readily enters solution and completely dissolves
- Easily penetrates organic soils
- Combination of anionic and nonionic detergents, which emulsifies blood and protein for quick removal
- Concentrated, economical use dilution ratio
- Free rinsing for film-free results
- Noncorrosive; safe on stainless steel, aluminum, gold and other soft metals
- Hard water tolerant, effective in all water conditions

Quality Statement: Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Usage Information

Water Hardness (Grains per gallon)	Product Dosage (oz./gal.)	(mL/L)
0 to 5	0.2	2.0
6 to 9	0.4	4.0
10 to 12	0.6	6.0
13+	1.0	8.0

Dilute product with warm water as shown in chart above. Soak instruments or utensils for 5 to 10 minutes. Drain and rinse. Process through instrument washer. Wear rubber gloves and safety glasses for prolonged contact.

If using a DEMA Proportioning Valve at the Sink — First Time Use with New Product

1. The metering tip may need to be replaced.
2. Acceptable colors are:
 - a. Tan (0.25 oz./gal.; 2 mL/L)
 - b. Orange (0.35 oz./gal.; 2.7 mL/L)
 - c. Turquoise (0.50 oz./gal.; 3.9 mL/L)
 - d. Pink (0.63 oz./gal.; 4.9 mL/L)
 - e. Clear (0.75 oz./gal.; 5.9 mL/L)
 - f. Brown (1 oz./gal.; 7.8 mL/L)

If using a Hand Pump at the Sink — As Needed

1. Each full pump equates to 1 oz. (29.6 mL).
 2. If specific sink volumes are unknown, the water line on the sink should be measured in order to determine proper dosing.
 3. See table above for sink dosage. For example, if the sink volume is six gallons, and the water hardness is ten grains per gallon, use three to four full pumps of product.
- **NOTE:** Discard solution if it becomes visibly soiled. Discard diluted solutions daily.
 - **NOTE:** Follow the instrument manufacturer's IFU for Decontamination/Cleaning/Sterilization.
 - **NOTE:** Follow hospital policy for appropriate PPE to be worn during reprocessing.
 - **NOTE:** Decontamination should be performed in accordance with AAMI ST79, AORN, and CDC guidelines. Please consult your Getinge Representative for additional information.
 - **CAUTION:** Do not add bleach, detergents, or any other chemical to this product. Added chemical may destroy the effectiveness of this cleaner.
 - **CAUTION:** Some plastics and rubbers might not be compatible with acid detergents or high water temperatures. Check with manufacturer for suitability.

Ordering Information

Getinge Clean	Order Number
3 x 5L	6036000000

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries.
Copyright 2021 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Intercon Chemical Co. 1100 Central Industrial Dr., St. Louis, MO 63110
Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA
Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada
www.getinge.com



MCV00044128 REV/B
EXP 11/23