



FLUOBEAM LM

Advanced technology to assess perfusion



FLUOBEAM LM

Advanced technology to assess perfusion

Over a decade of experience in fluorescence imaging, FLUOBEAM LM is an innovative solution for perfusion and lymphatic assessment.

FLUOBEAM LM with Indocyanine Green (ICG) is an integrated fluorescence imaging solution providing surgeons with a real-time image of the fluorescence in the operative field.

Its ease of use and ability to analyze images make it a major asset for surgeons.



Inject

Indocyanine Green (ICG)



Light

Lighting in situ of the fluorescent agent



Visualize

Real-time fluorescence assessment

FLUOBEAM LM

Panel PC (touch screen)
with its indication-oriented
software

Ergonomic optical head
Class 1 laser

Mobile cart
4 brakes

Control box

Limited footprint





Ergonomic and intuitive

The design of the FLUOBEAM LM optical head was developed to provide surgeons with a comfortable grip and intuitive handling during surgery.


Its joystick ensures simplified navigation in software functionalities and gives the surgeon total autonomy when using the device.

- Ease of use
- Intuitive navigation in the acquisition modes
- Ergonomic grip of the optical head
- A joystick to control software functionalities
- Complete autonomy of the surgeon from the sterile field

High performance and safety

FLUOBEAM LM has been designed to offer a high level of technical performance while ensuring the safety of users and patients.

The homogeneity of the laser excitation allows large areas to be imaged uniformly. This feature provides consistent and reliable information given by the displayed image.




Indication-oriented software

FLUOSOFT LM imaging software offers several acquisition modes depending on clinical indications, so that surgeons can get the expected images.

- Pseudo-colorization with push of a button to assess perfusion with relative quantification
- Automatic tool to help in perforator vessels identification before flap dissection
- Panoramic images of the lymphatic drainage

This approach, combined with proven acquisition protocols, helps prevent mishandling, and ensures great reliability.

FLUOBEAM LM in action



Plastic and Reconstructive surgery

Fluorescence imaging provides qualitative information and propose a relative quantification tool for better perfusion assessment. This information enables surgeons to carry out a real-time objective analysis of the perfusion quality.

During reconstructive surgery with autologous flaps (free or pedicle), perfusion assessment is essential to avoid postoperative complications such as total or partial skin or fat necrosis.

FLUOBEAM LM provides surgeons with accurate information to assess tissue perfusion during reconstructive surgery (breast, head and neck reconstruction procedures).

- Help in identifying intraoperatively perforator vessels and perforator angiosome
- Real-time intraoperative tissue perfusion assessment
- More precise flap design according to the perfused areas
- Relative quantification tool (additional information to improve the specificity of the method)
- Early identification of complications
- Postoperative monitoring

**Perforator
identification**



Free flap design



**Skin perfusion
assessment after
skin-sparing
mastectomy**



Axillary sentinel lymph node biopsy

The sentinel lymph node (SLN) procedure aims to detect, dissect and analyze lymph nodes. FLUOBEAM LM allows the visualization of superficial lymphatic drainage and the detection of sentinel lymph nodes after skin incision.

- High detection rate (similar to radioisotope)
- Direct visualization (lymphatic drainage and SLN)
- Simple patient care pathway (one injection only on the day of surgery)
- No exposure to radioactive product



Superficial lymphatic drainage visualization



Fluorescent sentinel lymph node in situ



Fluorescent sentinel lymph node ex vivo

Lymphatic drainage visualization

During consultation or during surgery, FLUOBEAM LM is used to image in real-time functional superficial lymphatic vessels. It helps surgeons to select the most appropriate treatment, especially for lymphedema patients.

- Early identification and localization of lymphatic drainage deficiencies to establish appropriate treatment
- Accurate visualization of functional lymphatic vessels before lymphaticovenous anastomosis (LVA)
- Accurate location of lymph nodes during lymph node transfer surgery
- Post-treatment monitoring
- Reverse mapping capability



Mapping of the lymphatic network



Visualization of small lymphatic vessels



Panoramic view of the whole limb



When two experts meet together

Keep your OR light on

Getinge Volista VisioNIR* allows to perform open surgeries using NIR fluorescence imaging in a seamless workflow by not turning on or off the surgical light.

Just keep the light on! No need to choose between one or the other.

- A powerful solution to guide surgeons, secure their actions with a better hand-eye coordination when using fluorescence guided surgery. No need to switch between on and off lighting.
- The surgical staff can stay focused on the patient on-going surgery. One less operation for the circulating staff.
- Uninterrupted workflow as you can keep the surgical light on during the entire procedure, no need to think about it.
- The ability to keep the OR light on provides better visibility of the operating room environment for the staff.
- Works simultaneously with the adjustable color temperature feature: while using Indocyanine Green (ICG) and NIR cameras, the surgeon can operate with the preferred color temperature. The dedicated enhancement mode improves the contrast on the screen.
- Keeping parameters of the OR Light like a standard mode with good color rendering, no change in shadow dilution or dimming.



One solution: unique filtered light

Thanks to the patented filters' wheel developed on Maquet Volista StandOP, the light emitted from the LEDs is filtered to reduce the remaining NIR wavelengths. Surgical lights disturbing the fluorescence signal emitted is now eliminated. Maquet Volista VisioNIR and NIR guided surgery cameras can be used simultaneously inside the operating room.



*available on Volista StandOP

More information on www.getinge.com



FLUOBEAM LM is a Class IIa medical device in Europe (CE 0197), manufactured by FLUOPTICS SAS.

FLUOBEAM LM is indicated to visualize on a screen the flow, the distribution and/or the accumulation of Indocyanine green (ICG) before, during and after surgery for the indications such as: visualization of the blood flow, visualization of the lymphatic flow, visualization and identification of the bile ducts during hepatobiliary surgery, visualization and detection of primary liver tumors and/or hepatic metastases.


This information is intended for an international audience outside the US.

This information is aimed exclusively at healthcare professionals or other professional audiences and is for informational purposes only, is not exhaustive and therefore should not be relied upon as a replacement of the Instructions for Use, service manual or medical advice. Getinge shall bear no responsibility or liability for any action or omission of any party based upon this material, and reliance is solely at the user's risk.

Any therapy, solution or product mentioned might not be available or allowed in your country. Information may not be copied or used, in whole or in part, without written permission by Getinge.

Views, opinions, and assertions expressed are strictly those of the interviewed and do not necessarily reflect or represent the views of Getinge.

Manufacturer · FLUOPTICS SAS · 44 rue des Berges, 38000 Grenoble, France · Phone: +33 4 85 87 06 71

© 2024 Getinge | Getinge and **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates | All rights reserved
VE01_FBLLM_06/24_EN_v2

www.getinge.com