





Intra-aortic Balloon Sizing Guide

 <p>< 5'0" (< 152 cm)</p>	 <p>5'0" to 5'4" (152 cm - 162 cm)</p>	 <p>5'4" to 6'0" (162 cm - 183 cm)</p>	 <p>5'4" & taller (≥ 162 cm)</p>
<p>30 cc</p>	<p>40 cc</p>	<p>50 cc</p>	
<p>Mega 7.5 Fr.</p>	<p>Sensation Plus 7.5 Fr. Mega 7.5 Fr.</p>	<p>Sensation Plus 8 Fr. Mega 8 Fr.</p>	
<p>25 cc</p>	<p>34 cc</p>	<p>40 cc</p>	
<p>Linear 7.5 Fr.</p>	<p>Sensation 7 Fr. Linear 7.5 Fr.</p>	<p>Sensation 7 Fr. Linear 7.5 Fr.</p>	

Note: This information is to be used as a guidance only. Clinical information and patient factors such as torso length should be considered when selecting the appropriate balloon size. Sensation and Sensation Plus are fiber-optic IAB catheters.

High Efficacy Intra-aortic Balloons

Why are the high efficacy IAB sizes better for my patients?


Larger volume balloons – Mega 30 cc, 40 cc, 50 cc and Sensation Plus 40 cc, 50 cc – displace more blood in the aorta during diastole, resulting in improved augmentation and unloading.

With these sizing guidelines, a larger patient population can benefit from **greater hemodynamic support** than ever before.

Benefits of larger volume IABs

- More blood volume displacement
- More diastolic augmentation
- More systolic unloading

Greater hemodynamic support and higher efficacy for every patient.



Hemodynamic Comparison*

* Bench testing completed by Getinge. Data on file. Bench test results are not necessarily predictive of clinical results.

PN: 0002-08-1092 Rev A · MCV00009481 REVD · *GETINGE , and Getinge, are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries · All rights reserved · CAUTION: Federal (U.S.A.) law restricts this device to sale, distribution and use by or on the order of a physician. ΔRefer to Instructions for Use for current indications, warnings, contraindications and precautions · 05/2018

Getinge • 1300 MacArthur Blvd., Mahwah, NJ 07430 · USA · +1 201 995-8700
45 Barbour Pond Drive, Wayne, NJ 07470 · USA · +1 888 880-2874

www.getinge.com

GETINGE 