Addendum to Servo Ventilator 300/300A - Operating Manual

Function test of control panel potentiometers and switches

The control panel potentiometers and switches must be tested at every 1000 and 3000 hour maintenance interval according to the following instructions.

Function test of control panel potentiometers

- Turn all potentiometer control knobs counter-clockwise to their end positions (the "Volume" control knob at least one turn).
- Set the mode selector to "Stand by".
- Initiate the Set Parameter Guide using the touchpad in a mode that will activate the potentiometer to be tested. A yellow light at the control knob shows if it is active.

Repeat the following procedure for all potentiometers:

- Turn the control knob clockwise to its end position (the "Volume" control knob at least one turn).
- Then turn the control knob slowly counter-clockwise to its end position (the "Volume" control knob at least one turn). While turning the control knob, check on the corresponding display or bargraph that:
 1. Displayed values reflect the position of the control knob.
 2. Displayed values follow the control knob setting evenly without sudden jumps.
 3. Displayed values or bargraph LEDs do not flash.

Function test of control panel switches

- Run the ventilator in "Pressure Control" mode with gases and a test lung connected.
- Turn the three switches "Oxygen breaths/Start breath", "Reset/2 min" and "Pause hold" back and forth to both positions. Check that the corresponding switch function is activated.
- Turn the three switches "Mode selector", "Patient range" and "Automode" step-by-step through all settings. Check that the yellow light at the corresponding setting is turned on.

If any malfunction is detected during these functions tests, the ventilator must not be connected to a patient. The malfunction must be remedied by MAQUET-authorized service personnel.

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