



Datasheet

**Servo-n**

**System version 4.1**

# Contents

|   |           |
|---|-----------|
| <b>Technical specifications .....</b>                 | <b>3</b>  |
| General .....   | 3         |
| The ventilator – general .....                        | 3         |
| Ventilation – general .....                           | 4         |
| User interface .....                                  | 4         |
| Power supply .....                                    | 4         |
| Gas supply .....                                      | 4         |
| Operating conditions .....                            | 5         |
| Non operating conditions .....                        | 5         |
| Standards – safety and functionality .....            | 5         |
| Communication / Interface .....                       | 5         |
| Invasive ventilation – modes .....                    | 6         |
| Invasive ventilation – leakage compensation .....     | 6         |
| Non invasive ventilation – modes .....                | 6         |
| Non invasive ventilation – leakage compensation ..... | 6         |
| High flow therapy (option) .....                      | 6         |
| VT/PBW .....  | 6         |
| Display .....   | 7         |
| Open Lung Tool trends (option) .....                  | 7         |
| Transpulmonary pressure measurement (option) .....    | 8         |
| Parameter settings .....                              | 8         |
| Backup parameter settings .....                       | 9         |
| Special functions .....                               | 9         |
| Disconnection .....                                   | 9         |
| Monitoring and trends .....                           | 9         |
| Alarms .....  | 10        |
| Autoset (alarm limits) specification .....            | 10        |
| Y sensor (option) .....                               | 11        |
| CO <sub>2</sub> analyzer (option) .....               | 11        |
| Edi module (option) .....                             | 12        |
| Aerogen nebulizers .....                              | 12        |
| Log function .....                                    | 12        |
| Saving of data .....                                  | 12        |
| Optional equipment .....                              | 13        |
| Compressor Mini (option) .....                        | 13        |
| Battery charger/calibrator (option) .....             | 13        |
| Service .....   | 13        |
| <b>Ordering information .....</b>                     | <b>13</b> |
| <b>Dimensional drawings .....</b>                     | <b>14</b> |
| Servo-n on mobile cart .....                          | 14        |
| Servo-n holder .....                                  | 14        |
| Servo-n (patient unit) on Servo-n holder .....        | 14        |

# Servo-n

## Technical specifications

### General

|                      |   |
|----------------------|---|
| Intended use         | The Servo-n ventilator system is: <ul style="list-style-type: none"> <li>intended for respiratory support, monitoring and treatment of neonatal and pediatric patients</li> <li>to be used only by healthcare providers</li> <li>to be used only in professional healthcare facilities and for transport within these facilities</li> </ul> |
| Clinical benefits    | Clinical benefits for Edi monitoring and NAVA: <ul style="list-style-type: none"> <li>to provide monitoring of the patient's breathing drive</li> <li>to improve synchrony between the ventilator system and patient when the electrical signal from the brain to the diaphragm is active</li> </ul>  |
| Instructions for use | Please carefully read the user's manual   |
| Legal manufacturer   | Maquet Critical Care AB   |
| Other products       | See separate data sheets.<br>Contact your local Getinge supplier for more information.  |

### The ventilator – general

|  | Servo-n   | Servo-n on mobile cart  |
|--|---|---|
| Base system weight                                   | Approximately 23 kg (50.7 lbs) <ul style="list-style-type: none"> <li>Patient unit 15 kg (33.0 lbs)</li> <li>User interface 4 kg (8.8 lbs)</li> <li>Handle 3 kg (6.6 lbs)</li> <li>Cable holder and cable 1 kg (2.2 lbs)</li> </ul> | Approximately 35 kg (77.2 lbs) <ul style="list-style-type: none"> <li>Base system approx. 23 kg (50.7 lbs)</li> <li>Mobile cart 13.0 kg (28.7 lbs)</li> </ul> |
| Dimensions of base (W x D), see dimensional drawings | 368 x 205 mm (14.5" x 8.1")   | 647 x 547 mm (25.5" x 21.5") incl. wheels   |
| Height (incl. user interface)                        | 826 mm (32.5")  | 1368 mm (53.8")   |
| Wheels   | N/A   | Four wheels with separate brakes  |
| A-weighted sound pressure level ( $L_{pA}$ )         | <40 dB, measured at a distance of 1 m (3.3 ft)  |   |
| A-weighted sound power level ( $L_{WA}$ )            | <51 dB  |   |

## Ventilation – general

|                              |   |
|------------------------------|---|
| Patient range                | <ul style="list-style-type: none"> <li>• Neonatal: 0.3–8 kg (0.7–17.6 lbs)</li> <li>• Pediatric: 3–30 kg (6.6–66.0 lbs)</li> </ul>  |
| Bias flow                    | 0.5 l/min *   |
| Internal compressible factor | Max. 0.1 ml/cmH <sub>2</sub> O  |
| Gas delivery system          | Microprocessor controlled valves  |
| Maximum airway pressure      | 125 cmH <sub>2</sub> O  |
| Method of triggering         | Flow, pressure and Edi (with Edi module and Edi catheter)   |
| Inspiratory flow range       | 0 to 33 l/min *   |
| Pressure drop                | <ul style="list-style-type: none"> <li>• Max. 6 cmH<sub>2</sub>O at a flow of 60 l/min (insp. channel)</li> <li>• Max. 3 cmH<sub>2</sub>O at a flow of 60 l/min (exp. channel)</li> </ul> |
| PEEP regulation              | Microprocessor controlled valve   |
| Expiratory flow range        | 0 to 192 l/min  |

\* In HFOV, adaptive BIAS flow and extended flow range.

## User interface

|              |  |
|--------------|--|
| Type         | TFT-LCD touchscreen                                    |
| Size         | 366 x 300 x 50 mm (14.4" x 11.8" x 2.0")               |
| Viewing area | 15" XGA, 1024 x 768 pixels with a 24-bit color palette |
| Weight       | Approximately 4 kg (8.8 lbs)                           |

## Power supply

|   |  |
|---|--|
| Power supply, automatic range selection   | <ul style="list-style-type: none"> <li>• 100–120 V AC, 2 A, 50–60 Hz</li> <li>• 220–240 V AC, 1 A, 50–60 Hz</li> </ul>   |
| Plug-in battery module:   |  |
| <ul style="list-style-type: none"> <li>• Battery backup (nickel-metal hydride, NiMH)</li> <li>• Battery capacity</li> <li>• Battery backup time</li> <li>• Recharge time</li> </ul> | <ul style="list-style-type: none"> <li>• Six battery module slots. Two batteries are delivered with the ventilator.</li> <li>• Rechargeable, 12 V, 3.5 Ah each</li> <li>• Ranging from 60 minutes (2 batteries) to 180 minutes (6 batteries)</li> <li>• Approximately 3 h/battery</li> </ul> |
| External 12 V DC  | 12.0 V–15.0 V DC, 10 A   |
| Typical min. power consumption (no optional modules, no ongoing battery charging, normal panel backlight)   | 100 VA, 40 W at 230 V or 75 VA, 40 W at 110 V  |
| Typical max. power consumption (with CO <sub>2</sub> , Edi and Y sensor modules, ongoing battery charging, max. panel backlight)  | 200 VA, 80 W at 230 V or 170 VA, 80 W at 110 V   |

## Gas supply

|  |   |
|--|---|
| Inlet gas pressure air/O <sub>2</sub>                    | 200–600 kPa / 2.0–6.0 bar / 29–87 PSI   |
| Connection standards available air/O <sub>2</sub>        | AGA, DISS, NIST, or French standard   |
| Inlet gas pressure HeO <sub>2</sub> (option)             | 340–600 kPa / 3.4–6.0 bar / 49–87 PSI   |
| Connection standards available HeO <sub>2</sub> (option) | AGA, DISS, NIST   |
| Unavailable gas/loss of gas pressure                     | The flow from an unavailable gas (air or O <sub>2</sub> ) is automatically compensated for so that the patient gets the preset volume and pressure. |
| Patient system gas connectors                            | Male 22 mm / female 15 mm. In accordance with ISO 5356-1  |
| Gas exhaust port   | Male 30 mm cone   |


## Operating conditions

|                                    |  |
|------------------------------------|--|
| Operating temperature              | +10 to +40°C (+50 to +104°F)   |
| Relative humidity                  | 15 to 95% non-condensing   |
| Atmospheric pressure               | 660 to 1060 hPa<br>(HFOV performance may be limited at high altitudes. See users instruction for details.) |
| Lowest pressure in patient circuit | -400 cmH <sub>2</sub> O  |

## Non operating conditions

|                                    |                              |
|------------------------------------|------------------------------|
| Temperature                        | -25 to +60°C (-13 to +140°F) |
| Relative humidity                  | <95% condensing              |
| Atmospheric pressure               | 470 to 1060 hPa              |
| Lowest pressure in patient circuit | N/A                          |

## Standards – safety and functionality

|  |   |
|--|---|
|                                 | The device complies with requirements and classification IIb of Medical Device Directive 93/42/EEC.<br><br>CE Mark Notified Body number: 0123.  |
| Classification   | IEC 60601-1: 2005 + A1:2012, Class I, continuous operation  |
| Standards  | <ul style="list-style-type: none"> <li>• ISO 80601-2-12:2011, ISO 80601-2-55:2018, EN 13544-1:2007+A1:2009</li> <li>• IEC 60601-1, Type B (equipment making physical contact with the patient and the gas pathways).</li> <li>• IEC 60601-1, Type BF (CO<sub>2</sub> analyzer, Y sensor, nebulizer patient unit and cable)</li> <li>• IEC 60601-1, Type CF defibrillation proof (Edi catheter and cable)</li> </ul> |
| Ingress protection   | IP 21   |
| Electromagnetic compatibility (EMC)  | According to limits specified in IEC 60601-1-2:2014   |
| The ' <b>Electromagnetic Compatibility Servo-u/Servo-n Ventilator System</b> ' is available from the manufacturer. |   |

## Communication / Interface

|                                     |   |
|-------------------------------------|---|
| Serial ports                        | <ul style="list-style-type: none"> <li>• Two RS-232C ports. For data communication via the Servo Communication Interface (SCI).</li> </ul>  |
| Servo Communication Interface (SCI) | A protocol for data communication with external devices   |
| Alarm output connection (option)    | <ul style="list-style-type: none"> <li>• 4-pin modular connector for communication of all active alarms</li> <li>• Switching capability: Max. 40 V DC, max. 500 mA, max. 20 W</li> </ul>  |
| Data transfer via USB port          | <ul style="list-style-type: none"> <li>• For transfer of trends, logs, screenshots and recordings to a USB memory stick</li> </ul>  |
| Ethernet port                       | <ul style="list-style-type: none"> <li>• The network connection (LAN) port is for service use, and should only be used by personnel trained and authorized by the manufacturer</li> </ul> |
| MSync, HL7 converter (optional)     | See separate datasheet  |
| VGA interface                       | VGA connector for duplication of the user interface   |

## Invasive ventilation – modes

|                        |  |
|------------------------|--|
| Controlled ventilation | <ul style="list-style-type: none"> <li>• PC (Pressure Control)</li> <li>• VC (Volume Control)</li> <li>• PRVC (Pressure Regulated Volume Control)</li> <li>• HFOV (Amplitude Control)</li> <li>• HFOV (V TGT) (Volume Target)</li> </ul> |
|------------------------|--|

Supported ventilation:

- PS/CPAP (Pressure Support / Continuous Positive Airway Pressure)
- VS (Volume Support)

|                   |  |
|-------------------|--|
| Automode (option) | <ul style="list-style-type: none"> <li>• Control mode: VC &lt;--&gt; Support mode: VS</li> <li>• Control mode: PC &lt;--&gt; Support mode: PS</li> <li>• Control mode: PRVC &lt;--&gt; Support mode: VS</li> </ul> |
|-------------------|--|

Combined ventilation

- SIMV (VC) + PS (Synchronized Intermittent Mandatory Ventilation)
- SIMV (PC) + PS
- SIMV (PRVC) + PS
- Bi-Vent/APRV (Airway Pressure Release Ventilation) (option)

|      |  |
|------|--|
| NAVA | <ul style="list-style-type: none"> <li>• Neurally Adjusted Ventilatory Assist via endotracheal tube or tracheostomy</li> </ul> |
|------|--|

VC and SIMV (VC) + PS and Automode VC <--> VS are not available in the neonatal patient category.

HFOV and HFOV (V TGT) are only available in neonatal patient category

## Invasive ventilation – leakage compensation

|                                 |  |
|---------------------------------|--|
| Max. leakage compensation level | <ul style="list-style-type: none"> <li>• Neonatal: - 25 l/min</li> </ul> |
|---------------------------------|--|

## Non invasive ventilation – modes

|                        |  |
|------------------------|--|
| Controlled ventilation | <ul style="list-style-type: none"> <li>• NIV PC (option)</li> </ul>  |
| Supported ventilation  | <ul style="list-style-type: none"> <li>• NIV PS (option)</li> <li>• Nasal CPAP</li> </ul>                                    |
| NIV NAVA               | <ul style="list-style-type: none"> <li>• Neurally Adjusted Ventilatory Assist via non-invasive patient interfaces</li> </ul> |

NIV PS is not available in the neonatal patient category

## Non invasive ventilation – leakage compensation

|                                 |  |
|---------------------------------|--|
| Max. leakage compensation level | <ul style="list-style-type: none"> <li>• Pediatric and neonatal: <ul style="list-style-type: none"> <li>- Inspiratory: up to 33 l/min</li> <li>- Expiratory: up to 25 l/min</li> <li>- Nasal CPAP: up to 20 l/min</li> </ul> </li> </ul> |
|---------------------------------|--|

Disconnection flow (configurable)

- Low:
  - 7.5 l/min
- High:
  - 15 l/min
- Disabled: the ventilator system will continue to deliver assist even when leakage is excessive.

|                      |                                   |
|----------------------|-----------------------------------|
| Connection detection | Manual or automatic via bias flow |
|----------------------|-----------------------------------|

## High flow therapy (option)

|                    |   |
|--------------------|---|
| Flow setting range | <ul style="list-style-type: none"> <li>• Pediatric: 0.5–30 l/min</li> <li>• Neonatal: 0.5–20 l/min</li> </ul> |
|--------------------|---|

## VT/PBW

|                             |  |
|-----------------------------|--|
| Predicted Body Weight (PBW) | Automatically calculated for adult patients based on gender and height (130–200 cm)                              |
| Body Weight (BW)            | Entered for neonatal and pediatric patients, as well as adult patients shorter than 130 cm or taller than 200 cm |
| VT/PBW (VT/BW) in ml/kg     | Automatically calculated, displayed and trended  |

## Display

|                     |  |
|---------------------|--|
| Views               | <ul style="list-style-type: none"> <li>• Basic</li> <li>• Advanced</li> <li>• Loops</li> <li>• Servo Compass (option)</li> <li>• Pes &amp; PL (option)</li> <li>• Distance</li> <li>• Family</li> </ul> <p>Each of the screen layout views offers a specific combination of displayed waveforms, loops and presented values.</p> |
| Real time waveforms | <ul style="list-style-type: none"> <li>• Airway Pressure</li> <li>• Flow</li> <li>• Volume</li> <li>• Edi (option)</li> <li>• CO<sub>2</sub> (option)</li> <li>• Transpulmonary pressure (option)</li> <li>• Esophageal pressure (option)</li> </ul>   |
| Loops               | <ul style="list-style-type: none"> <li>• Pressure – Volume</li> <li>• Pressure – Flow</li> <li>• Volume – Flow</li> </ul> <p>A reference loop and three overlaying loops can be displayed.</p>   |
| Servo Compass       | Visualizes volume (VT/BW) and pressure (total or driving) in relation to set targets in invasive modes.  |
| Short trends        | <ul style="list-style-type: none"> <li>• During ventilation in all ventilation modes, short trends of the numerical values in the first column can be displayed.</li> <li>• Trend time 15 minutes to 72 hours.</li> </ul>  |
| Trends              | <ul style="list-style-type: none"> <li>• Trending of measured and calculated values.</li> <li>• Trend time 1 to 72 hours.</li> <li>• Order of trended values can be set by the user.</li> </ul>  |

## Open Lung Tool trends (option)

### OLT trends (option)

|                       |   |
|-----------------------|---|
| Graphical trend areas | <p>1:</p> <ul style="list-style-type: none"> <li>- Pei (end-inspiratory pressure)</li> <li>- Pdrive *</li> <li>- PEEP</li> </ul> <p>2:</p> <ul style="list-style-type: none"> <li>- VT<sub>CO<sub>2</sub></sub> (when applicable)</li> <li>- SI * (Stress Index, adult patient category only)</li> <li>- C<sub>dyn</sub></li> </ul> <p>3 (standard):</p> <ul style="list-style-type: none"> <li>- VT<sub>i</sub></li> <li>- VT<sub>e</sub></li> </ul> <p>3 (option):</p> <ul style="list-style-type: none"> <li>- PL<sub>ei</sub></li> <li>- PL<sub>ee</sub></li> <li>- PL<sub>drive</sub> *</li> </ul> |
|-----------------------|---|

\* Pdrive, PL drive and SI only shown as values – not graphical trends

|                       |  |
|-----------------------|--|
| Modes                 | All invasive modes   |
| Trend time            | 5, 10, 15, 30, or 60 minutes                                       |
| Recruitment recording | Recording of recruitments for retrospective review of recruitments |

### Auto RM (option)

|   |   |
|---|---|
| Automatic recruitment maneuver with two phases for adult and pediatric patients | Available in PC, PRVC and VC invasive ventilation modes |
|---|---|

|                 |  |
|-----------------|--|
| Maneuver phases | <p>1. Recruitment made in PC mode with I:E set to 1:1. PEEP and inspiratory pressure increase according to a preset pattern.</p> <p>2. Post-recruitment, where the system returns to the mode set prior to recruitment and sets a user-selected post-recruitment PEEP.</p> |
|-----------------|--|

|                        |  |
|------------------------|--|
| Recruitment parameters | <ul style="list-style-type: none"> <li>• PEEP<sub>max</sub></li> <li>• RR</li> <li>• P<sub>max</sub></li> <li>• Δ PEEP/step</li> <li>• Breaths/step</li> <li>• Breaths at P<sub>max</sub></li> <li>• Post-RM PEEP</li> </ul> |
|------------------------|--|

|                      |  |
|----------------------|--|
| Recruitment analysis | Pre- and Post-recruitment measurements during 5 breaths each |
|----------------------|--|

|                       |   |
|-----------------------|---|
| Recruitment recording | Automatic recording of recruitments with retrospective review of recruitments possible in OLT trends or as recruitment recordings |
|-----------------------|---|

## Transpulmonary pressure measurement (option)

|   |   |  |
|---|---|--|
| Esophageal pressure measurement via Auxiliary pressure (Paux) port on Y sensor module |   |  |
| Pes Catheter Positioning  | Automatic maneuver to validate Esophageal balloon positioning and filling |  |
| Waveforms   | Pes   | Esophageal pressure  |
|   | PL  | Transpulmonary pressure = Paw – Pes  |
| Numerical values  | PL ei   | End inspiratory PL = Paw ei – Pes ei   |
|   | PL ee   | End expiratory PL = PEEP – Pes ee  |
|   | PL drive  | PL ei – PL ee (passive ventilation)<br>PL max (inspiration) – PL ee (active breathing)                               |
|   | ΔPes  | Pes max (inspiration) – Pes ee (positive Pes deflection)<br>Pes min (inspiration) – Pes ee (negative Pes deflection) |

## Parameter settings

| Parameter   | Neonatal range | Pediatric range |
|---|----------------|-----------------|
| Tidal volume (ml)                                       | 2–50           | 10–350          |
| Minute volume (l/min)                                   | 0.1–7.5        | 0.3–20          |
| Apnea, time to alarm (s)                                | 1–45           | 2–45            |
| Max. apnea time in Automode (s)                         | 3–15           | 3–15            |
| Pressure level above PEEP (cmH <sub>2</sub> O)          | 0–79           | 0–79            |
| Pressure level above PEEP (in NIV) (cmH <sub>2</sub> O) | 0–60           | 0–60            |
| P <sub>mean</sub> (HFOV) (cmH <sub>2</sub> O)           | 5–40           | –               |
| P <sub>amp</sub> (HFOV) (cmH <sub>2</sub> O)            | 0–100          | –               |
| Tidal Volume (HFOV) (ml)                                | 0.2–40         | –               |
| I:E <sub>ratio</sub> (HFOV)                             | 1:3–1:1        | –               |
| Frequency (HFOV) (Hz)                                   | 5–20           | –               |
| PEEP (cmH <sub>2</sub> O)                               | 0–50           | 0–50            |
| PEEP in NIV (cmH <sub>2</sub> O)                        | 2–20           | 2–20            |
| CPAP pressure (cmH <sub>2</sub> O)                      | 2–20           | 2–20            |
| Respiratory rate (breaths/min)                          | 4–150          | 4–150           |
| SIMV rate (breaths/min)                                 | 1–60           | 1–60            |
| Breath cycle time, SIMV (s)                             | 0.5–15         | 0.5–15          |
| P <sub>High</sub> (cmH <sub>2</sub> O)                  | 2–50           | 2–50            |
| T <sub>High</sub> (s)                                   | 0.2–30         | 0.2–30          |
| T <sub>PEEP</sub> (s)                                   | 0.1–10         | 0.1–10          |
| PS above P <sub>High</sub> (cmH <sub>2</sub> O)         | 0–78           | 0–78            |
| O <sub>2</sub> concentration (%)                        | 21–100         | 21–100          |
| I:E ratio   | 1:10–4:1       | 1:10–4:1        |
| Ti (s)  | 0.1–5          | 0.1–5           |
| NAVA level (cmH <sub>2</sub> O/μV)                      | 0–15           | 0–15            |
| Edi trigger (μV)  | 0.1–2.0        | 0.1–2.0         |
| T <sub>Pause</sub> (s)                                  | –              | 0–1.5           |
| T <sub>Pause</sub> (% of breath cycle time)             | –              | 0–30            |
| Flow trigger (l/min)                                    | 0–0.5          | 0–0.5           |
| Pressure trigger (cmH <sub>2</sub> O)                   | -1 to -20      | -1 to -20       |
| Insp. rise time (% of breath cycle time)                | 0–20           | 0–20            |
| Insp. rise time (s)                                     | 0–0.2          | 0–0.2           |
| End inspiration (% of peak flow)                        | 1–70           | 1–70            |
| End inspiration (% of peak flow) in NIV                 | 10–70          | 10–70           |
| Decelerating flow pattern in VC (%)                     | --             | 0–100           |
| Flow adaptation in VC                                   | --             | on/off          |



## Backup parameter settings

| Parameter  | Neonatal range | Pediatric range |
|--|----------------|-----------------|
| Inspiratory tidal volume (ml)                                | 2–50           | 10–350          |
| Pressure level above PEEP in backup (cmH <sub>2</sub> O)     | 5–79           | 5–79            |
| Pressure level above PEEP in NIV backup (cmH <sub>2</sub> O) | 5–60           | 5–60            |
| Respiratory rate in backup (breaths/min)                     | 4–150          | 4–150           |
| I:E ratio  | 1:10–4:1       | 1:10–4:1        |
| Ti (s)   | 0.1–5          | 0.1–5           |

## Special functions

| Special function              | Setting range  |
|-------------------------------|--|
| Manual breath                 | Initiation of 1 breath (In SIMV mode initiation of 1 mandatory breath)                           |
| Static measurements           | Insp. or exp. hold (0–30 seconds)  |
| Nebulization                  | 5–30 min/Continuous/Off  |
| O <sub>2</sub> boost level    | Off, 1–79 %, 100%  |
| O <sub>2</sub> boost function | Activate O <sub>2</sub> boost up to 1 minute   |
| Leakage compensation          | On/Off   |
| Circuit compensation          | On/Off   |
| Edi monitoring                | In all ventilation modes, in High Flow therapy and in Standby (with Edi module and Edi catheter) |
| Previous mode                 | Activates previously used mode   |
| Backup ventilation            | Backup On/Off  |
| Apnea management              | Several parameters   |
| Pause oscillation             | Pause oscillation during HFOV  |

## Disconnection

|                         |   |
|-------------------------|---|
| Pre-oxygenation time    | Max. 2 min                                |
| Post-oxygenation time   | Max. 1 min                                |
| Patient disconnected    | High priority alarm activated after 1 min |
| Adjustable oxygen level | 21 – 100 %                                |

## Monitoring and trends

|  |                                       |
|--|---------------------------------------|
| Peak airway pressure                         | Ppeak                                 |
| Pause airway pressure                        | Pplat                                 |
| Mean airway pressure                         | Pmean                                 |
| Driving airway pressure                      | Pdrive                                |
| Positive end expiratory pressure             | PEEP                                  |
| Continuous positive airway pressure          | CPAP                                  |
| Spontaneous breaths per minute               | RR sp                                 |
| Respiratory rate                             | RR                                    |
| Spontaneous expiratory minute volume         | MVe sp                                |
| Inspired minute volume                       | MVi                                   |
| Expired minute volume                        | MVe                                   |
| Leakage fraction (%)                         | Leakage                               |
| Inspired tidal volume                        | VTi                                   |
| Expired tidal volume                         | VTe                                   |
| End expiratory flow                          | Flowee                                |
| Measured oxygen concentration                | O <sub>2</sub> conc                   |
| CO <sub>2</sub> end tidal concentration      | etCO <sub>2</sub>                     |
| CO <sub>2</sub> minute elimination           | VCO <sub>2</sub>                      |
| CO <sub>2</sub> tidal elimination            | VTCO <sub>2</sub>                     |
| CO <sub>2</sub> Diffusion (HFOV)             | DCO <sub>2</sub>                      |
| Dynamic compliance                           | Cdyn                                  |
| Static compliance                            | Cstatic                               |
| Inspiratory resistance                       | Ri                                    |
| Expiratory resistance                        | Re                                    |
| Work of breathing, ventilator                | WOBvent                               |
| Work of breathing, patient                   | WOBpat                                |
| Elastance                                    | E                                     |
| P 0.1  | P 0.1                                 |
| Shallow Breathing Index                      | SBI                                   |
| Peak Edi value                               | Edipeak                               |
| Average Edipeak                              | Edipeak average (monitoring only)     |
| Average Edimin                               | Edimin average (monitoring only)      |
| Minimum Edi value                            | Edimin                                |
| Ratio of expired tidal volume to body weight | VT/BW                                 |
| Switches to backup per minute                | Backup $\Sigma$ (trended value only)  |
| Time in backup in percent per minute         | Backup % (trended value only)         |
| Tidal Volume (HFOV)                          | VT <sub>HF</sub>                      |
| Pressure Amplitude (HFOV)                    | P <sub>amp</sub>                      |
| I:E-ratio (HFOV)                             | I:E <sub>HF</sub>                     |
| Heliox gas consumption                       | HeO <sub>2</sub> (trended value only) |

## Alarms

| Alarm  | Neonatal range  | Pediatric range              |
|--|---|------------------------------|
| Airway pressure (upper alarm limit)                    | 16–90 cmH <sub>2</sub> O  | 16–90 cmH <sub>2</sub> O     |
| Airway pressure NIV (upper alarm limit)                | 16–70 cmH <sub>2</sub> O  | 16–70 cmH <sub>2</sub> O     |
| Respiratory rate (upper alarm limit)                   | 2–160 breaths/min   | 2–160 breaths/min            |
| Respiratory rate (lower alarm limit)                   | 1–159 breaths/min   | 1–159 breaths/min            |
| Expired minute volume (upper alarm limit)              | 0.02–30 l/min   | 0.02–30 l/min                |
| Expired minute volume (lower alarm limit)              | 0.01 – 20 l/min   | 0.01 – 20 l/min              |
| End expiratory pressure (upper alarm limit)            | 1–55 cmH <sub>2</sub> O   | 1–55 cmH <sub>2</sub> O      |
| End expiratory pressure (lower alarm limit)            | Off, 1–47 cmH <sub>2</sub> O  | Off, 1–47 cmH <sub>2</sub> O |
| No patient effort (Apnea) alarm                        | 1–45 s  | 2–45 s                       |
| Automatic return to support mode on patient triggering |   |                              |
| No consistent patient effort                           | Yes, described in User's manual   |                              |
| High continuous pressure                               | Yes, described in User's manual   |                              |
| O <sub>2</sub> concentration                           | Set value ±5 vol% or ≤18 vol%   |                              |
| Gas supply   | Below 200 kPa (2.0 bar/29 PSI), above 600 kPa (6.0 bar/87 PSI)  |                              |
| Battery  | <ul style="list-style-type: none"> <li>Limited battery capacity: 10 min.</li> <li>No battery capacity: less than 3 min</li> <li>Low battery voltage.</li> </ul> |                              |
| End tidal CO <sub>2</sub> (upper and lower limit)      | 0.5–19.9 %, 4–99mmHg, 0.5–13.9kPa   |                              |
| Leakage too high                                       | Yes, described in User's manual   |                              |
| Technical  | Yes, described in User's manual   |                              |
| Mean Airway Pressure (HFOV) High / Low                 | Yes   | –                            |
| Pressure amplitude (HFOV) High / Low                   | Yes   | –                            |
| VT <sub>HF</sub> High / Low                            | Yes   | –                            |
| VT <sub>HF</sub> Limited                               | Yes   | –                            |

## Autoset (alarm limits) specification

| Autoset (alarm limits) specification                        | Invasive ventilation, controlled modes only                                 |
|---|---|
| High airway pressure  | Mean peak pressure +10 cmH <sub>2</sub> O or at least 35 cmH <sub>2</sub> O |
| Inspiratory tidal volume too high                           | The greater of VTi + 30 % or VTi +2 ml                                      |
| Expiratory minute volume (upper alarm limit)                | Mean expiratory minute volume +50 %   |
| Expiratory minute volume (lower alarm limit)                | Mean expiratory minute volume -50 %   |
| Respiratory rate (upper alarm limit)                        | Mean respiratory rate +40 %   |
| Respiratory rate (lower alarm limit)                        | Mean respiratory rate -40 %   |
| End expiratory pressure (upper alarm limit)                 | Mean end expiratory pressure +5 cmH <sub>2</sub> O                          |
| End expiratory pressure (lower alarm limit)                 | Mean end expiratory pressure -3 cmH <sub>2</sub> O                          |
| End tidal CO <sub>2</sub> concentration (upper alarm limit) | Mean end tidal CO <sub>2</sub> concentration +25 %                          |
| End tidal CO <sub>2</sub> concentration (lower alarm limit) | Mean end tidal CO <sub>2</sub> concentration -25 %                          |
| High frequency Tidal Volume in HFOV (upper alarm limit)     | High frequency Tidal Volume +22 % *   |
| High frequency Tidal Volume in HFOV (lower alarm limit)     | High frequency Tidal Volume -29 % *   |
| Pressure amplitude in HFOV (V TGT) (upper alarm limit)      | Mean Pressure amplitude in HFOV +40 % *                                     |
| Pressure amplitude in HFOV (V TGT) (lower alarm limit)      | Mean Pressure amplitude in HFOV -40 % *                                     |

\* These alarm limits correspond to a 50 % increase / decrease of DCO<sub>2</sub>

## Y sensor (option)

| Y sensor (option)       | Size  | Weight          |
|-------------------------|---|-----------------|
| Y sensor module         | W 154 x L 90 x H 21 mm<br>(W 6.1" x L 3.5" x H 0.8")  | 280 g (0.6 lbs) |
| Y sensor                | W 18 x L 50 x H 27 mm<br>(W 0.7" x L 2.0" x H 1.1")   | 11 g            |
| Connectors and cables   | <ul style="list-style-type: none"> <li>• 15 mm male and female conical connector on flow sensor according to ISO 5356-1</li> <li>• Pressure port on module, pressure line, 2.0 m (6.6 ft), phthalate free PVC</li> <li>• Flow sensor cable, 2.0 m (6.6 ft)</li> </ul> |                 |
| Sensor material         | • Single use: PC, Polycarbonate   |                 |
| Power source            | Powered by the ventilator system, $\leq 4.5$ W during normal operation  |                 |
| Measuring method        | Hot Wire Anemometer (HWA)   |                 |
| Parameters              | <ul style="list-style-type: none"> <li>• Airway pressure</li> <li>• Airway flow</li> <li>• Inspiratory and expiratory volumes</li> <li>• Trigger and End inspiration</li> </ul>   |                 |
| Measuring range         | <ul style="list-style-type: none"> <li>• Flow: 0.12 to 32 l/min</li> <li>• Pressure: -40 to 120 cmH<sub>2</sub>O</li> </ul>   |                 |
| Y sensor resistance     | 10 cmH <sub>2</sub> O/l/s at 30 l/min   |                 |
| Dead space              | $\leq 1$ ml   |                 |
| Pressure line connector | Gable mounted bulk head connector to fit tubing with an inner diameter of 3-4 mm (0.12–0.16")   |                 |

Y sensor is recommended for conventional ventilation with small tidal volumes, and required during HFOV to get flow and tidal volumes measurements.

## CO<sub>2</sub> analyzer (option)

| CO <sub>2</sub> analyzer (option)    | Size   | Weight               |
|--------------------------------------|--|----------------------|
| CO <sub>2</sub> analyzer module      | W 154 x L 90 x H 21 mm<br>(W 6.1" x L 3.5" x H 0.8")   | 265 g<br>(0.58 lbs)  |
| Sensor (Capnostat 5)                 | 32.0 x 47.0 x 21.6 mm<br>(1.3" x 1.9" x 0.8")  | 20 g                 |
| Operating temperature                | 10 to 33 °C (50 to 91 °F)  |                      |
| Airway adapter                       | 10 g   |                      |
| Power source                         | Powered by the ventilator  |                      |
| Connectors and cables                | Sensor   | 2.8 m (9.2 ft) cable |
| Measuring method                     | Mainstream, dual-wavelength, non-dispersive infrared   |                      |
| Parameters                           | <ul style="list-style-type: none"> <li>• CO<sub>2</sub> end tidal concentration (etCO<sub>2</sub>)</li> <li>• CO<sub>2</sub> minute elimination (VCO<sub>2</sub>)</li> <li>• CO<sub>2</sub> tidal elimination (VTCO<sub>2</sub>)</li> </ul>    |                      |
| Measuring range                      | <ul style="list-style-type: none"> <li>• 0 to 100 mmHg CO<sub>2</sub> partial pressure</li> <li>• 0 to 13.3 kPa CO<sub>2</sub> partial pressure</li> <li>• 0 to 13.2 % CO<sub>2</sub> volume (at a barometric pressure of 1013 hPa)</li> </ul> |                      |
| System response time CO <sub>2</sub> | The total system response time of the CO <sub>2</sub> monitor when exposed first to air and then to a gas mix with 5.0 % CO <sub>2</sub> is <250 ms  |                      |
| Warm-up time                         | 15 s to initial CO <sub>2</sub> indication maximum 2 minutes to full specification   |                      |
| Oxygen concentration compensation    | Automatic. Values supplied from the ventilator system  |                      |
| Barometric pressure compensation     | Automatic. Values supplied from the ventilator system  |                      |
| Digitizing rate                      | 100 Hz   |                      |
| Airway adapter dead space            | • Neonatal/pediatric: <1 cm <sup>3</sup>   |                      |

## Edi module (option)

| Edi module (option) | Size   | Weight               |
|---------------------|--|----------------------|
| Edi module          | W 154 x L 90 x H 21 mm<br>(W 6.1" x L 3.5" x H 0.8")   | 0.25 kg<br>(0.6 lbs) |
| Edi catheter cable  | 2.0 m (6.6 ft)   | -                    |
| Power source        | Powered by the ventilator  |                      |
| Power consumption   | <3 W during normal operation   |                      |
| Parameters          | <ul style="list-style-type: none"> <li>• Edi waveform</li> <li>• ECG leads waveforms</li> <li>• NAVA estimated pressure waveform (Pedi)</li> </ul> |                      |

## Log function

|                |   |
|----------------|---|
| Event log      | <ul style="list-style-type: none"> <li>• Alarms</li> <li>• Ventilator settings</li> <li>• Apnea periods</li> <li>• Maneuvers and O<sub>2</sub> boost</li> </ul>                               |
| Diagnostic log | <ul style="list-style-type: none"> <li>• Technical alarms</li> <li>• Test results</li> <li>• Service records</li> <li>• Software installation</li> <li>• Configuration information</li> </ul> |

## Aerogen nebulizers

| Aerogen nebulizers | Pro   | Solo  |
|--------------------|---|---|
| Size               | W 50 x L 50 x H 45 mm<br>(W 2.0" x L 2.0" x H 1.8") | W 48 x L 25 x H 67 mm<br>(W 1.9" x L 1.0" x H 2.6") |
| Weight             | Approx. 25 g (0.88 oz)                              | Approx. 14 g (0.49 oz)                              |
| Particle size      | 1–5 µm mass median aerodynamic diameter (MMAD)      |   |
| Flow rate          | >0.2 (average: ~0.4) ml/min                         |   |
| Max. volume        | 10 ml   | 6 ml  |
| Residual volume    | <0.1 ml for 3 ml dose                               |   |
| Control cable      | 1.8 m (5.9 ft)                                      |   |

## Saving of data

|  |   |
|--|---|
| Recording of current waveform and parameter values | 30 seconds of data will be recorded (15 seconds before and 15 seconds after activation). Up to 40 recordings can be stored. |
| Saving screenshots                                 | Up to 40 screenshots can be stored.   |
| Saving recruitments                                | Up to 12 manual and/or automatic recruitment recordings can be stored (option).   |
| Export files                                       | Recordings, screenshots, recruitments, trends and event log can be saved and exported to a USB memory stick.                |

## Optional equipment

| Optional equipment          | Weight                | Dimensions  | Maximum load              |
|-----------------------------|-----------------------|---|---------------------------|
| Mobile cart                 | 15.0 kg<br>(33.1 lbs) | W 647 x L 547 x H 557 mm<br>(W 25.5" x L 21.5" x H 21.9")   | -                         |
| Drawer for mobile cart      | 0.6 kg<br>(1.3 lbs)   | W 247 x L 118 x H 302 mm<br>(W 9.7" x L 4.6" x H 11.9")   | -                         |
| Shelf base                  | 2.5 kg<br>(5.5 lbs)   | W 207 x L 302 x H 43 mm<br>(W 8.2" x L 4.6" x H 1.7")   | -                         |
| Pendant/bed holder          | 3.2 kg<br>(7.1 lbs)   | W 302 x L 302 x H 393 mm<br>(W 11.9" x L 11.9" x H 15.5")   | -                         |
| Humidifier holder           | 0.6 kg<br>(1.3 lbs)   | W 243 x L 38 x H 185 mm<br>(W 9.6" x L 1.5" x H 7.3")   | 5 kg<br>(11.0 lbs)        |
| Support arm 178             | 2.2 kg<br>(4.8 lbs)   | L 900 mm (35.4")  | 1–3 kg<br>(2.2–6.6 lbs) * |
| * depending on angle        |                       |   |                           |
| User interface holder       | 0.6 kg<br>(1.3 lbs)   | W 46 x L 90 x H 123 mm<br>(W 1.8" x L 3.5" x H 4.8")  | -                         |
| Cable holder for handle     | 0.1 kg<br>(0.2 lbs)   | W 138 x L 92 x H 155 mm<br>(W 5.4" x L 3.6" x H 6.1")   | 5 kg<br>(11.0 lbs)        |
| Waterbag/IV pole            | 0.4 kg<br>(0.9 lbs)   | W 148 x L 26 x H 1007 mm<br>(W 5.8" x L 1.0" x H 39.6")   | 1.5 kg<br>(3.3 lbs)       |
| Gas cylinder restrainer kit | 1.0 kg<br>(2.2 lbs)   | Upper:<br>W 104 x L 65 x H 48 mm<br>(W 4.1" x L 2.5" x H 1.9")<br>Lower:<br>W 106 x L 162 x H 76 mm<br>(W 4.1" x L 6.4" x H 3.0") | Two 4.5 liter bottles     |
| Y piece holder              |                       | W 26 x L 52 x H 46 mm<br>(W 1.0" x L 2.0" x H 1.8")   |                           |

## Compressor Mini (option)

See separate datasheet

## Battery charger/calibrator (option)

See separate datasheet

## Service


|                     |  |
|---------------------|--|
| Regular maintenance | Preventive maintenance must be performed by authorized personnel at least once every 5000 hours of operation or once every 12 months, whichever comes first. |
|---------------------|--|

## Ordering information


Servo-n, ventilator system and accessories:  
See separate information in "System flowchart, Servo-n"

# Dimensional drawings


## Servo-n on mobile cart



## Servo-n holder



## Servo-n (patient unit) on Servo-n holder







Getinge is a global provider of innovative solutions for operating rooms, intensive care units, sterilization departments and for life science companies and institutions. Based on our firsthand experience and close partnerships with clinical experts, healthcare professionals and medtech specialists, we are improving the everyday life for people, today and tomorrow.

Servo-n may be pending regulatory approvals to be marketed in your country. Contact your Getinge representative for more information. This document is intended to provide information to an international audience outside of the US.

May be changed without prior notice.

**Manufacturer** · Maquet Critical Care AB · Röntgenvägen 2 SE-171 54 Solna · Sweden · +46 (0)10 335 73 00

**[www.getinge.com](http://www.getinge.com)**