

XEN

Chemistries
for the ultimate
peace of mind

Instrument Lubricant+



Cart Washer



Manual



Washer/
Disinfector

Instructions for Use

XEN Instrument Lubricant+ is a water-based liquid lubricant for use in automated washer applications, including washer-disinfectors and cart washers, as well as manual use. Instrument Lubricant+ lubricates the joints and moving parts of medical instruments, reducing friction corrosion and contributing to a longer life for the equipment. It does not interfere with steam or ethylene oxide sterilization.

Order no.	Description
6036004701	Instrument Lubricant+ 1 x 5 L.
6036004801	Instrument Lubricant+ 1 x 5 gal.
6036004901	Instrument Lubricant+ 1 x 15 gal.



Directions for Use in Automated Washing Equipment

1. Instrument Lubricant+ should be automatically dosed into the final rinse cycle. Dosing is 0.05–0.13 oz./gal. (0.4–1.0 ml/L). Automated washers should be set to dose at >122°F (>50°C). Deionized (DI) or reverse osmosis (RO) water is recommended for the final rinse.



Directions for Use Manual Method

1. In a suitable container, mix 20-40ml/L in clean potable water, preferably deionized RO or distilled at <77°F (25°C).
2. Open instrument joints and place into basket.
3. Submerge instruments in the solution for 30 seconds.
4. Remove instruments from solution and allow to drain and dry prior to sterilization.

The solution must be discarded daily or after each shift, or if it becomes visibly contaminated.

GETINGE 

Caution

- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.

Danger: Causes serious eye damage.

First Aid: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a poison center or doctor. For specific treatment refer to section 4 of SDS.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Contains polyalkylene glycol, benzoisothiazolinone. For more ingredient information, visit www.getinge.com

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Instrument Lubricant+	•	•	•	•	•	•	•	•	•	•			○	○	•	○

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instruments/device manufacturer

Teflon™ is a trademark of Chemours.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kingdom.
Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA
Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

