

Instrument Lubricant+



Cart Washer



Manual



Washer/
Disinfectant

Product Specification

Description

A water-based, silicone- and mineral-oil-free lubricant for use the final rinse phase in manual or automated washer applications. Instrument Lubricant+ lubricates joints and moving parts of medical instruments, reducing friction corrosion and prolonging instrument life. It will not interfere with steam or ethylene oxide sterilization processes.

Application

Instrument Lubricant+ is used in the final rinse phase to lubricate and protect instruments. Refer to equipment and instrument/device manufacturers' instructions before use.

Key Features

- Silicone- and mineral-oil-free—Steam and EtO permeable, so does not interfere with sterilization or affect the functionality of general surgical instruments
- Clear solution—Translucency increases technician safety with sharps
- Safe for staff—Does not carry the aspiration toxicity classification and has been tested according to ISO 10993 for biocompatibility and found not to be irritating, sensitizing or cytotoxic
- Reduces friction and wear—May aid in prolonging the lifespan of surgical instruments and help reduce the downtime and costs associated with friction, stress/corrosion, and maintenance/repair
- Environmentally friendly—Biodegradable formulation and recyclable packaging

Ingredients

Contains nonionic surfactants, anionic surfactants <5 %, and preservative.



Technical Data

Specific Gravity/Density (68°F/20°C)	1.004 (typical)
pH as supplied (typical)	9.0
pH at 0.26 oz./gal. (2 ml/L) water (typical)	8.6
Viscosity (concentrate, 68°F/20°C)	4.2 cSt (typical)
Expiration/Use By Date	36 months from date of manufacture in original sealed containers
Dosing (Manual)	2.56–5.12 oz./gal. (20–40 ml/L)
Dosing (Automated)	0.05–0.13 oz./gal. (0.4–1.0 ml/L)
Product Dimensions	5 L = 7.6" L x 5.7" W x 10" H (193 x 146 x 255 mm)
	5 gal. = 11.6" L x 9.4" W x 15.4" H (295 x 240 x 390 mm)
	15 gal. = 15.7" Dia. x 24" H (399 x 609 mm)

Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome-plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Instrument Lubricant+	•	•	•	•	•	•	•	•	•	•			○	○	•	○

• = The product is suitable for the material

○ = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Ordering Information

Size	Order Number
1 x 5 L	6036004701
1 x 5 gal.	6036004801
1 x 15 gal.	6036004901

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE 

MCV00110900 REVA