

XEN

Chemistries
for the ultimate
peace of mind

Instrument Shine



Cart Washer



Maintenance



Washer/
Disinfectant

Product Specification

Description

An all-in-one cleaner for brightening stainless steel, removing stains and mineral scale, and protecting it from corrosion by repassivating the instrument's surface to prolong its useful lifetime. Use in soaking applications or automated washers, including washer-disinfectors and cart washers.

Application

For revitalizing of stainless steel medical instruments and automated washers. Refer to equipment and instrument/device manufacturer's instructions before use.

Key Features

- Recommended for cleaning stainless steel instruments
- Improves instrument performance and prolongs life-span—Reduces corrosion and removes contaminants from all parts of stainless steel instruments
- Removes rust stains, discoloration, and hard water scale from stainless steel instruments, equipment, and automated washer chambers—Brightens, destains, and improves performance of equipment
- Restores stainless steel finish—Repassivates instruments, prolonging lifespan
- Environmentally friendly—Biodegradable formulation and recyclable packaging

Ingredients

Orthophosphoric acid and nitric acid.



Technical Data

| | |
|--|--|
| Specific Gravity/Density (68°F/20°C) | 1.28 (typical) |
| pH as supplied (typical) | <1 |
| pH at 0.26 oz./gal. (2 ml/L) water (typical) | 2 |
| Viscosity (concentrate, 68°F/20°C) | 5.4 cSt (typical) |
| Expiration/Use By Date | 36 months from date of manufacture in original sealed containers |
| Temperature Range (Automated) | 160°F (71°C) |
| Dosing (Automated) | 0.5–1.0 quart for every 8 gal. |
| Dosing (Manual) | 2 fl oz. per quart of hot water (50ml/L) |
| Product Dimensions | 1 qt. = 34.8"Dia x 8.8"H (885 x 223 mm) |

GETINGE 

Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

| | Aluminum | Anodized Aluminum | Titanium | Copper | Brass | Gold | Nickel | Cobalt Chrome | Matte chrome-plated | Stainless steel | Glass | Ceramics | Rubber | Plastic | Teflon™/ PTFE | Polycarbonate |
|------------------|----------|-------------------|----------|--------|-------|------|--------|---------------|---------------------|-----------------|-------|----------|--------|---------|---------------|---------------|
| Instrument Shine | | | | | | | | | | • | • | | | | | |

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

Storage Information

Store in a cool, dry place protected from frost and away from alkalis and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

Cautions

- Do not use on polycarbonate material
- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PP.
- Refer to Safety Data Sheet (SDS) for further safety and handling

Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Ordering Information

| Size | Order Number |
|-----------|--------------|
| 6 x 1 qt. | 6036004501 |

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kindom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE

MCV00110902 REVA