

# XEN

Chemistries  
for the ultimate  
peace of mind

## Neutral



Ophthalmic  
Instruments



Orthopedics



Robotic  
Instruments



Flexible  
Scopes



Rigid  
Scopes



Cart Washer



Manual



Ultrasonic



Washer/  
Disinfectant

## Product Specification

### Description

A pH-neutral and nonenzymatic detergent, making it well-suited for ophthalmic instrumentation. It is low foaming to optimize performance of and prevent damage to automated washers, ultrasonics, and lumen and robotic instruments.

### Application

For cleaning of robotic and ophthalmic instruments, endoscopes, and other devices by manual or automated methods. Refer to equipment and instrument/device manufacturers' instructions and AAMI ST79 and AORN cleaning recommendations before use.

### Key Features

- Suitable for manual cleaning, automated cleaning, and ultrasonics
- pH-neutral—Recommended by AAMI ST79 and AORN to be compatible with the widest range of instrument surfaces
- Low foaming—Optimizes performance of and prevents damage to automated washers and ultrasonic cleaners; prevents damage to delicate lumen or hollow instruments or devices
- Nonenzymatic—Safe and gentle for ophthalmic instruments; eliminates the risk of enzymes causing Toxic Anterior Segment Syndrome (TASS)
- Efficient surfactant system—Good cleaning performance without the risk of enzyme-related complications
- Exceptional material compatibility—Allows multiple types of instruments to be processed together, increasing throughput

### Ingredients

Surfactants, corrosion inhibitors, and preservatives



### Technical Data

Specific Gravity/Density (68°F/20°C)	1.038 (typical)
pH as supplied (typical)	8.8
pH at 0.64 oz./gal. (5 ml/L) water (typical)	8.6
Viscosity (concentrate, 68°F/20°C)	4.5 cSt (typical)
Expiration/Use By Date	24 months from date of manufacture in original sealed containers
Temperature Range (Automated)	112°F–140°F (50°C–60°C)
Dosing (Automated)	0.13–2.05 oz./gal. (1.0–16.0 ml/L)
Temperature Range (Manual and Sonic)	77°F–113°F (25°C–45°C)
Dosing (Manual and Sonic)	0.13–2.05 oz./gal. (1.0–16.0 ml/L)
Product Dimensions	5 L = 7.6"L x 5.7"W x 10"H (193 x 146 x 255 mm)
	5 gal. = 11.6"L x 9.4"W x 15.4"H (295 x 240 x 390 mm)
	15 gal. = 15.7"Dia. x 24"H (399 x 609 mm)

GETINGE 

## Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte-chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Neutral	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

## Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

## Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

## Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

## Ordering Information

Size	Order Number
1 x 5 L	6036003201
1 x 5 gal.	6036005701
1 x 15 gal.	6036005801

## Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kindom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

[www.getinge.com](http://www.getinge.com)

**GETINGE** 

MCV00110897 REVA