

XEN

Chemistries
for the ultimate
peace of mind

Sanizyme+



Orthopedics



Robotic
Instruments



Flexible
Scopes



Rigid
Scopes



Ultrasonic

Product Specification

Description

A recharged formulation of the well-known original, Sanizyme+ is a low-foaming, pH-neutral, triple-enzyme detergent optimized for the cleaning of medical instruments and devices in ultrasonics. It is effective against a wide range of organic contaminants and compatible with a wide range of materials.

Application

For cleaning of surgical-grade instruments/reusable items, endoscopes, and other devices. Refer to equipment and instrument/device manufacturers' instructions and AAMI ST79 and AORN cleaning recommendations before use.

Key Features

- Suitable for ultrasonic washers—The neutral pH prevents damage to the interior chamber, while low foam prevents interference with the ultrasonics
- pH-neutral—Recommended by AAMI and AORN to prevent pitting and staining of medical instruments
- Low foaming—Reduces chance of damage to the instruments or interference with the ultrasonics
- Triple-enzyme formulation—Protease enzymes break down proteins in blood and tissue; lipase enzymes attack fats and lipids; and amylase enzymes break down carbohydrates
- Good in adverse water conditions—Efficiently penetrates, surrounds, and carries away surface contamination
- Exceptional material compatibility—Allows a wide variety of instruments in the same load

Ingredients

Protease, lipase, and amylase enzymes, surfactants, sequestering agents, and corrosion inhibitors.



Technical Data

Specific Gravity / Density (68°F/20°C)	1.039 (typical)
pH as supplied (typical)	8.7
pH at 0.64 oz./gal. (5 ml/L) water (typical)	8.4
Viscosity (concentrate, 68°F/20 °C)	3.8 cSt (typical)
Expiration/Use By Date	24 months from date of manufacture in original sealed containers
Temperature Range	77°F–113°F (25°C–45°C)
Dosing	1.02 oz./gal. (8 ml/L)
Product Dimensions	5 L = 7.6"L x 5.7"W x 10"H (193 x 146 x 255 mm)
	5 gal. = 11.6"L x 9.4"W x 15.4"H (295 x 240 x 390 mm)
	15 gal. = 15.7"Dia. x 24"H (399 x 609 mm)

GETINGE 

Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte-chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Sanizyme+	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

• = The product is suitable for the material

◦ = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Ordering Information

Size	Order Number
1 x 5 L	6036005201
1 x 5 gal.	6036005301
1 x 15 gal.	6036005401

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andresseway Way, Derby, DE21 6YH, United Kingdom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE

MCV00110904 REVA