

XEN

Chemistries
for the ultimate
peace of mind


Universal



Orthopedics



Robotic
Instruments



Rigid
Scopes



Cart Washer



Washer/
Disinfector

Instructions for Use

XEN Universal is a mild alkaline liquid detergent for use in automated washer applications, including washer-disinfectors and cart washers. Universal is low foaming to optimize performance of the automated washer. It is safe to use on robotic instruments. Its three enzymes—protease, lipase, and amylase—break down blood and proteins, fats and oils (lipids), and carbohydrates, respectively.

Order no.	Description
6036003601	Universal 1 x 5 L.
6036003701	Universal 1 x 5 gal.
6036003801	Universal 1 x 15 gal.



Directions for Use in Automated Washing Equipment

1. The detergent should be automatically dosed into the wash cycle. Dosing is 0.13–0.51 oz./gal. (1–4 ml/l), depending on soil load and water quality. Higher doses will be required for higher soil loads and for poor water quality. Automated washers should be set to dose at 95°F (35°C). The wash temperature should be set at 122°F–140°F (50°C–60°C).

GETINGE 

Caution

- Do not mix with other chemicals.
- Wear correct PPE and follow applicable guidelines. Additional information on SDS.
- Check material compatibility before use.

Danger: May cause allergy or asthma symptoms or breathing difficulties if inhaled. Avoid breathing fumes, mists, vapors, or spray. In case of inadequate ventilation, wear respiratory protection.

First Aid: IF INHALED: If breathing is difficult, remove person to fresh air and keep at rest in a comfortable position for breathing. If experiencing respiratory symptoms, call a poison center or doctor.

Dispose of contents and container in accordance with local, regional, and international regulations.

Contains triethanolamine, protease, amylase, and lipase enzymes. For more ingredient information, visit www.getinge.com

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Universal	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instruments/device manufacturer

Teflon™ is a trademark of Chemours.

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kingdom.
Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA
Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada
www.getinge.com

