

XEN

Chemistries
for the ultimate
peace of mind

Xcelerate+



Ophthalmic Instruments



Orthopedics



Robotic Instruments



Flexible Scopes



Rigid Scopes



Pretreatment

Product Specification

Description

A ready-to-use, nonenzymatic, pH-neutral pretreatment foam spray for point-of-use and reprocessing delays to prevent drying of surgical soils and enhance the cleaning process. The wide spray angle makes it simple to get complete coverage. Suitable for stainless steel, ceramics, glass, plastic, and soft metal surfaces. Water-based, it is easily rinsed off.

Application

For point-of-use (prior to transport) on medical instruments and devices, including ophthalmic equipment. Begin the decontamination process as soon as possible to prevent drying of surgical soils, post procedure per AAMI ST79 and AORN recommendations. Refer to instrument/device manufacturer's instructions before use.

Key Features

- Ready to Use—Quick and easy to use. Product does not require mixing, eliminating exposure
- Automatic spray—Quickly covers all instruments in the tray
- pH-neutral—Recommended by AAMI/AORN for cleaning delicate medical instruments and devices; reduces chance of pitting or staining
- Nonenzymatic—Safe for use on ophthalmic instruments and outperformed enzymatic sprays in studies
- Keeps soils moist—Prevents drying and hardening for up to 100 hours in open-air conditions
- Foam spray—Allows easy verification of instrument coverage and breaks down within 10 minutes to improve visibility and avoid sharps injury
- Free-rinsing—Rinses cleanly for greater residue reduction



Ingredients

Surfactants, rheology modifier, preservative and corrosion inhibitors.

Technical Data

Specific Gravity/Density (68°F/20°C)	1 (typical)
pH as supplied (typical)	8
Viscosity (concentrate, 68°F/20°C)	6 cSt (typical)
Expiration/Use By Date	24 months from date of manufacture in original sealed containers
Temperature Range	Designed for use at ambient temperature
Product Dimensions	17 oz. = 25.6" Dia. x 9.8"H (650 x 250 mm)

GETINGE

Material Compatibility

Always follow the cleaning instructions from the manufacturer of the items and instruments and check material for compatibility before use. The following chart serves only as a guide.

	Aluminum	Anodized Aluminum	Titanium	Copper	Brass	Gold	Nickel	Cobalt Chrome	Matte chrome-plated	Stainless steel	Glass	Ceramics	Rubber	Plastic	Teflon™/ PTFE	Polycarbonate
Xcelerate+	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

- = The product is suitable for the material
- = Pay special attention to the guidelines from the instrument/device manufacturer

Teflon™ is a trademark of Chemours.

Storage Information

Store in a cool, dry place protected from frost and away from acids and strong oxidizing agents. Store upright in original containers. Optimum storage temperature 41°F–77°F/5°C–25°C

Cautions

- Check material for compatibility before use
- Do not mix with other chemicals
- Always wear correct PPE
- Refer to Safety Data Sheet (SDS) for further safety and handling

Environmental Awareness

Biodegradable formulation and recyclable packaging.

Dispose of contents and container in accordance with local, regional, national, and international regulations.

Ordering Information

Size	Order Number
12 x 17 oz.	6036005101

Service

For technical support in regards to appropriate detergent use, please contact your local Getinge representative.

Or click or scan here:



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. Copyright 2023 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kingdom

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE 

MCV00110890 REVA