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Références

1. Vignaux L, Grazioli S, Piquilloud L, Bochaton N, Karam O, Jaecklin T, Levy-jamet Y, Tourneux P, Jolliet P, Rimensberger P. Optimizing patient ventilator synchrony during invasive ventilator assist in children and infants remains a difficult task. *Ped Crit Care Med*. 2013;14(7): 316-325.
2. Beck J, Reilly M, Grasselli G, et al. Patient-ventilator interaction during neurally adjusted ventilatory assist in low birth weight infants. *Pediatr Res* 2009;65(6):663-8.
3. Longhini F, Ferrero F, De Luca D, et al. Neurally adjusted ventilatory assist in preterm neonates with acute respiratory failure. *Neonatology*. 2015;107(1):60-7
4. Mally PV, Beck J, Sinderby C, et al. Neural breathing pattern and patient-ventilator interaction during neurally adjusted ventilatory assist and conventional ventilation in newborns. *Pediatr Crit Care Med* 2018;19(1):48-55.
5. Beck J, Emeriaud G, Liu Y, Sinderby C. Neurally-adjusted ventilatory assist (NAVA) in children: a systematic review. *Minerva Anesthesiol* 2016;82(8):874-83.
6. Beck J, Sinderby C. Neurally adjusted ventilatory assist in newborns. *Clin Perinatol* 2021, Dec;48(4):783-811.
7. Tabacaru CR, Moores Jr RR, Khoury J, Rozycki HJ. NAVA-synchronized compared to nonsynchronized noninvasive ventilation for apnea, bradycardia, and desaturation events in VLBW infants. *Pediatr Pulmonol*. 2019 Nov;54(11):1742-6
8. Hovespyan K, Firestone KS, Moore J, Stein H. Effects of NAVA Compared to SIMV Ventilation on Cardiac Function in Preterm Neonates. *Resp Care* 2020;65(10):3451491.
9. Surkov D M. Neurally Adjusted Ventilatory Assist Mode of Mechanical Ventilation in Neonates with Hypoxic-Ischemic Encephalopathy. *Ukr J Cardio Surg* 2019;2(35):102-11.
10. Baudin F, Emeriaud G, Essouri S, Beck J, Javouhey E, Guerin C. Neurally adjusted ventilatory assist decreases work of breathing during non-invasive ventilation in infants with severe bronchiolitis. *Crit Care* 2019;23(1):120.
11. De la Oliva P, Schuffelmann C, Gomez-Zamora A, Vilar J, Kacmarek RM. Asynchrony, neural drive, ventilatory variability and COMFORT: NAVA vs pressure support in pediatric patients. A randomized cross-over trial. *Int Care Med*. 2012;38:838-46.
12. Piastra M, De Luca D, Costa R, et al. Neurally adjusted ventilatory assist vs pressure support ventilation in infants recovering from severe acute respiratory distress syndrome: nested study. *J Crit Care* 2014;29(2):312:e1-5.
13. Firestone KS, Beck J, Stein H. Neurally Adjusted Ventilatory Assist for Noninvasive Support in Neonates. *Clin Perinatol*. 2016 Dec;43(4):707-24.
14. Chidini G, De Luca D, Calderini E, et al. Implementation of noninvasive neurally adjusted ventilatory assist in pediatric acute respiratory failure: a controlled before-after quality improvement study. *J Anesth Analg & Crit Care* 2021(Sep);1:1-9
15. Kallio M, Peltoniemi O, Anttila E, Pokka T, Kontiokari T. Neurally Adjusted Ventilatory Assist (NAVA) in Pediatric Intensive Care – A Randomized Controlled Trial. *Pediatr Pulmonol*. 2015 Jan; 50(1):55-62.
16. Lee J, Kim HS, Jung YH, et al. Neurally adjusted ventilatory assist for infants under prolonged ventilation. *Pediatr Int* 2017;59(5):540-4
17. Baez Hernandez N, Milad A, Li Y, et al. Utilization of neurally adjusted ventilatory assist (NAVA) mode in infants and children undergoing congenital heart surgery: a retrospective review. *Pediatr Cardiol* 2019;40(3):563-9.
18. Gibu CK, Cheng PY, Ward RJ, et al. Feasibility and physiological effects of noninvasive neurally adjusted ventilatory assist in preterm infants. *Pediatr Res* 2017; 82(4):650-7.
19. Parikka V, Beck J, Zhai Q, Leppäsalo J, Lehtonen L, Soukka H. The effect of caffeine citrate on neural breathing pattern in preterm infants. *Early Hum Dev*. 2015 Oct;91(10):565-8.
20. Rahmani A, Ur Rehman N, Chedid F. Neurally adjusted ventilatory assist (NAVA) mode as an adjunct diagnostic tool in congenital central hypoventilation syndrome. *J Coll Physicians Surg Pak* 2013; Feb;23(2):154-156
21. Makker K, Cortez J, Jha K, et al. Comparison of extubation success using noninvasive positive pressure ventilation (NIPPV) versus noninvasive neurally adjusted ventilatory assist (NI-NAVA). *J Perinatol* 2020;40(8):1202-10.
22. Keszler, M. Volume-targeted ventilation: one size does not fit all. Evidence based recommendations for successful use. *Arch Dis Child Fetal Neonatal Ed* 2018;1-5
23. Lee BK, Shin SH, Jung YH, et al. Comparison of NIV-NAVA and NCPAP in facilitating extubation for very preterm infants. *BMC Pediatr* 2019 Aug 28;19(1):298
24. Morita et al. The usability of ventilators: a comparative evaluation of use safety and user experience *Crit Care* (2016) 20:263
25. Ducharme-Crevier L, Du Pont-Thibodeau G, Emeriaud G. Interest of Monitoring Diaphragmatic Electrical Activity in the Pediatric Intensive Care Unit. *Crit Care Res Pract*. 2013; 384210.
26. Stein H, Firestone K. Application of neurally adjusted ventilatory assist in neonates. *Semin Fetal Neonatal Med*. 2014 Feb;19(1):60-9.
27. Emeriaud G, Larouche A, Ducharme-Crevier L et al. Evolution of inspiratory diaphragm activity in children over the course of the PICU stay. *Intensive Care Med*. 2014 Nov;40(11):1718-26
28. Protain AP, Firestone KS, McIninch NL, Stein HM. Evaluating peak inspiratory pressures and tidal volume in premature neonates on NAVA ventilation. *Eur J Pediatr*. 2021 Jan;180(1):167-175
29. Parikka V, Beck J, Zhai Q, Leppäsalo J, Lehtonen L, Soukka H. The effect of caffeine citrate on neural breathing pattern in preterm infants. *Early Hum Dev*. 2015 Oct;91(10):565-8
30. Amigoni A, Rizzi G, Divisic A et al. Effects of propofol on diaphragmatic electrical activity in mechanically ventilated pediatric patients. *Intensive Care Med*. 2015 Oct;41(10):1860-1.
31. Lee, J, Parikka V, Lehtonen, L, Soukka H. Parent-infant skin-to-skin contact reduces the electrical activity of the diaphragm and stabilizes respiratory function in preterm infants. 2021 *Pediatr Res* June.
32. Baudin F, Emeriaud G, Essouri S et al. Physiological effect of prone position in children with severe bronchiolitis: a randomized cross-over study (BRONCHIO- DV). *J Pediatr* 2019;205:112-9.
33. Wolf G, Walsh B, Green M, Arnold J. Electrical activity of the diaphragm during extubation readiness testing in critically ill children. *Pediatr Crit Care Med* 2010;12:e220e4.
34. Poets CF, Roberts RS, Schmidt B, et al. Canadian Oxygen Trial Investigators. Association between Intermittent Hypoxemia or Bradycardia and Late Death or Disability in Extremely Preterm Infants. *JAMA*. 2015 Aug 11;314(6):595-603
35. Morley CJ, Davis PG, Doyle LW et al. COIN Trial Investigators. Nasal CPAP or intubation at birth for very preterm infants. *N Engl J Med*. 2008 Feb 14;358(7):700-8.
36. SUPPORT Study Group of the Eunice Kennedy Shriver NICHD Neonatal Research Network, Finer NN, et al. Early CPAP versus surfactant in extremely preterm infants. *N Engl J Med*. 2010 May 27;362(21):1970-9.
37. Courtney SE, Pyon KH, Saslow JG, et al. Lung recruitment and breathing pattern during variable versus continuous flow nasal continuous positive airway pressure in premature infants: an evaluation of three devices. *Pediatrics*. 2001 Feb;107(2):304-8.
38. Dargaville PA, Gerber A, Johansson S, et al. Incidence and Outcome of CPAP Failure in Preterm Infants. *Pediatrics*. 2016;138:e20153985-e20153985.
39. Houtekie L, Moerman D, Bourleau A, et al. Feasibility Study on Neurally Adjusted Ventilatory Assist in Noninvasive Ventilation After Cardiac Surgery in Infants. *Respir Care*. 2015 Jul;60(7):1007-14
40. Longhini F, Scarlino S, Gallina MR et al. Comparison of neurally-adjusted ventilator assist in infants before and after extubation. *Minerva Pediatr* 2018 Apr;70(2):133-40
41. Yagui ACZ, Goncalves PA, Murakami SH, et al. Is noninvasive neurally adjusted ventilatory assistance (NIV-NAVA) an alternative to NCPAP in preventing extubation failure in preterm infants? *J Matern Fetal Neonatal Med* 2019;1-151
42. Shetty S, Evans, K, Cornuauud P, et al. Neurally Adjusted Ventilatory Assist in very prematurely born infants with evolving/established bronchopulmonary dysplasia. *AJP Rep* 2021;11(4):e127-e131.
43. Così G, Monzani A, Genoni G et al. Neurally Adjusted Ventilatory Assist (NAVA): a prospective interventional study in neonates. *Minerva Pediatr*. 2019 Apr 05
44. Stein H, Howard D. Neurally Adjusted Ventilatory Assist (NAVA) in Neonates less than 1500 grams: a retrospective analysis. *J Pediatr* 2012;160:786e9.
45. Meyers, N, Rodrigues N, Ari A. High-frequency oscillatory ventilation: A narrative review. *Can J Res Ther* 2019. 55:40-6
46. Bordessoule A, Piquilloud L, Lyazidi A, Moreira A, Rimensberger PC. Imposed Work of Breathing During High-Frequency Oscillatory Ventilation in Spontaneously Breathing Neonatal and Pediatric Models. *Resp Care* 2018 Sep. 63(9):1085-1093.
47. Galindo-Filho, V.C. et al. Radioaerosol Pulmonary Deposition Using Mesh and Jet Nebulizers During Noninvasive Ventilation in Healthy Subjects. *Respir. Care* 2015, 60(9):1238-1246

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