

Data Notice – Cardiosave

Purpose of this notice

The purpose of this data notice is to inform you, in accordance with the transparency requirements for connected products under Article 3.2 of the EU Data Act, about the product data generated through your use of Cardiosave and your rights regarding such data.

Data that is accessible to you through the device's user interface and typically utilized as part of providing care is not included in this data notice. Also note that the information below regarding type of data, format, volume etc. may vary based on the use and setup in each individual case. Please contact us if you have any questions related to your devices.

Product data generated by the connected product

During use, the device may generate the following types of data:

Type: Log data: Recent Triggers, Recent Alarms, Electrical Test Failures, Recent Faults, Power Activity Log, and Battery Statistics.

Format: Hexadecimal data

Estimated volume: All data (currently 980 bytes) is sent in a single packet every 500ms.

Capable of generating data continuously and in real time: Yes

Storage of product data and retention time

The product data will be stored in Cardiosave (on device in Non-Volatile Memory (NVRAM) and on a remote server, if such connecting has been set up by the user.

Retention time for product data in Cardiosave depends on the usage rate, since capacity is limited. Data is retained until it is overwritten.

You can delete product data in Cardiosave at any time.

How you may access the product data

Product data can be accessed directly by the user in the non-clinical service mode through exporting it to a USB drive.

Contact details for access requests

To ask questions about data access, contact Getinge through:

E-mail: data.protection@getinge.com