

Accessories

For Getinge Thoracic Chest Drainage Systems



19910 — Red Cap

- For capping patient line before disposal
- For use if only one patient line is needed on a dual drain
- Twenty-five (25) each per case — sterile



19913 — Adult/Pediatric Connector 1/4 inch x 3/8 inch

- For use with an adult drain and smaller catheters — 16, 20 French
- Twenty-five (25) each per case — sterile



19914 — Standard Adult Connector 3/8 inch x 3/8 inch

- Standard 5-in-1 connector
- For use on larger catheters — 24 French and larger
- Twenty-five (25) each per case — sterile



19915 — Pediatric Connector 1/4 inch x 1/4 inch

- Larger connector included with infant/pediatric drain
- For use with an adult drain and smaller catheters
 - 16 French or larger thoracic catheters, 9 and 11 mm mediastinal catheters
- Twenty-five (25) each per case — sterile



19916 — Sims Connector 3/16 inch x 1/4 inch


- Same connector included with infant/pediatric drain
- For use with smaller catheter — 8, 10, 12 French
- Twenty-five (25) each per case — sterile



19918 — Step Down Connector 3/16 inch x 3/8 inch

- For use with adult chest drain
- Fits smaller catheters — 8, 10, 12 French
- Twenty-five (25) each per case — sterile

Accessories: For Getinge Thoracic Chest Drainage Systems



19907 — Patient Tube Pack

- 3/8 inch ID tubing
 - Used if patient tube requires changing; i.e. sterility is accidentally compromised
 - Ten (10) each per case — sterile
-



19928 — Suction Y-Adapter

- For use when connecting 2 chest drains to 1 suction source
 - Ten (10) each per case — sterile
-



19930 – Express™ Suction Bulb

- For use with the Express™ Mini 500 (No. 16400)
 - Allows portable, manual vacuum application
 - Ten (10) each per case — sterile
-

None of our drain accessories are made with natural rubber latex

Manufactured by Atrium Medical Corporation / 40 Continental Blvd., Merrimack, NH 03054 / Tel. (603) 880-1433 · Protected by the following international and U.S. patent(s): <http://patents.maquet.com>. · Getinge, GETINGE, Atrium and Express are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries · Getinge and Atrium are registered with the U.S. Patent and Trademark Office. · Copyright 2022 Atrium Medical Corp. · All rights not expressly granted are reserved. · Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. · PN 010112 Rev AF 12/22

www.getinge.com

This document is intended to provide information to an international audience outside of the US.

GETINGE 