

Green Tape

CI 159/.75, CI 159/1

SteriTec Green Tape tape adhesive is specifically formulated to adhere to disposable wrap. After the steam sterilization process, the green indicator lines will turn black/dark brown indicating that the wrapped unit was processed.



Key Features

- 4 minutes at 132°C and 3 minutes at 135°C
- Not made with natural rubber latex or lead
- Process indicator turns black when exposed to proper steam sterilization conditions

Performance Characteristics

Indicator stripes change from light green to black

Storage

The storage conditions are printed on the product packaging

Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

FDA 510(K)

K161024

Standards Classification

- FDA Sterilization Process Indicator
- ANSI/AAMI ST79, Type 1
- ISO 11140-1 Type 1

SteriTec Part No.

Packaging

CI 159/.75, Green Tape 3/4 inch	48 rolls
CI 159/1, Green Tape 1 inch	36 rolls