

EO Sterilization Indicator Tape

CI 124/.75, CI 124/1

The EO Sterilization Indicator Tape is for use in ethylene oxide sterilizers operating at 600 mg/l EO gas, 40-60% relative humidity, 130°F, for 45 minutes or longer. The word “gas” turns green after exposure to a EO gas sterilization process, thus providing identification of processed items.



Key Features

- Process indicator turns from brown to green when exposed to proper ethylene oxide sterilization conditions

Performance Characteristics

Tape must adhere to item to be sterilized throughout the sterilization process. The process indicator word “gas” on the tape will change color, brown to green when exposed to ethylene oxide gas.

Color Change

Brown to green

Storage

The storage conditions are printed on the product packaging. Store away from oxidizing agents

Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

FDA 510(K)

K002861

Standards Classification

- FDA Process Indicator
- ISO 11140-1, Type 1

SteriTec Part No.

CI 124/.75,
EO Sterilization Tape 3/4"

Packaging

48 rolls

CI 124/1,
EO Sterilization Tape 1"

36 rolls