


Container Card

CI 138

The Container Card chemical indicators are designed for use in steam sterilizers operating at 121°C or 134°C. When used as directed, the Container Card indicator gives visible confirmation that the container was exposed to steam.



Key Features

- 121°C or 134°C
- Can be used as a permanent record

Performance Characteristics

When exposed to steam sterilization conditions, the bar changes color. The card can fit in various container card-holders by folding or not folding the top sizing flap to the appropriate size.

Physical Characteristics

Dimensions: L 70mm x W 39mm

Color Change

White to black

Stated Values

- 10 minutes at 121°C
- 2 minutes at 134°C

Storage

The storage conditions are printed on the product packaging

Shelf Life

36 months from date of manufacture. The expiry date is printed on the product packaging

Standards Classification

ISO 11140-1, Type 1

SteriTec Part No.

CI 138, Container Card

Packaging

8 boxes of 500