

Getinge Clean

Detergent Viscosities

Name	Viscosity		Density lb/gal
	cP	cS	
PC Acid Det.	16	24	12.6
PC Alkaline Det.	10	14	11.3
PC Rinse Add.	20	21	8.6
PC Neutral	6	7	9.99
PC AI Safe	12	16	10.95
PC 3x Enzyme	4	4	8.35
Neutrawash	2	2	8.9
Neutrawash Plus	3	3	8.9
Acid Detergent	3	3	9.08
Alkaline Detergent	4	5	11
Acid Neutralizer	6	7	9.88
Tec Cage Wash	5	6	10.2
Glass Cage & Rack	2	2	8.8
Tec Rinse	20	21	8.6
Tec Wash III	4	5	9.7
Instrument Lube	8	8	8.55
Instrument Lube Plus	4	4	8.16
Renuzyme	2	2	9.1
Renuzyme Plus	5	5	8.76
Renuzyme WR	2	2	9.1

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries.
Copyright 2020 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Sales Office, US • Getinge • 1 Geoffrey Way • Wayne, NJ 07470

Sales Office, Canada • 90 Matheson Blvd West, Suite 300 • Mississauga, Ontario L5R 3R3 • Canada

www.getinge.com

MCV00106251 REVA
EXP07/22