

# EZTest EO Biological Indicators

## EZEO

EZTest® (Ethylene Oxide) is for biological testing of ethylene oxide sterilizers.

The Self-Contained Biological Indicator EZTest contains *Bacillus Atrophaeus*.



### Key Features

- *Bacillus Atrophaeus*
- 48 hour results with no laboratory required

### Performance Characteristics

After exposure to ethylene oxide, the biological indicator, when incubated, will indicate either growth or no growth.

### Color Change

- BI Vial Process Indicator - Blue to brown
- BI Growth Medium - Red to yellow with growth

### Physical Characteristics

- Dimensions: L 48mm x W 10mm

### Incubation Temperature/Time

35°-39°C, 48 hours

### Storage

The storage conditions are printed on the product packaging. Do not store near sterilants or other chemicals

### Shelf Life

24 months from date of manufacture. The expiry date is printed on the product packaging

### 510(K)

Mesa Labs, K930683

### Standards Classification

Biological Indicator - ANSI/AAMI/ISO 11138

### SteriTec Part No.

EZEO, EZTest EO Biological Indicator

### Packaging

Box of 100