

Getinge Clean Instrument Lubricant

Ensures lubrication of hinges, pins and joints

Low dosage, bio-compatible lubricant that extends instrument life

Getinge Clean Instrument Lubricant is added in the final rinse phase for instrument lubrication. Its low dosage emulsion formulation ensures deposition of a thin film of lubricant.

When used regularly, Getinge Clean Instrument Lubricant can improve instrument life reducing the long-term financial impact of premature instrument replacement.

This bio-compatible lubricant does not harden during sterilization and does not inhibit the sterilization process.



Available in 5L, 10L and 25L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at getinge.com/contact