

Getinge Clean Pre-treatment Foam

— Ready-to-use pre-treatment foam spray

Applied on medical devices immediately after use to prevent drying of surgical soils

The Getinge Clean Pre-treatment Foam is available in 500 ml hand-held dispensers. It is used principally for spraying surgical instruments and other utensils after use, before being transferred to the CSSD. It is especially recommended when there is a delay between instrument use and decontamination.

Pre-treatment Foam is not intended for machine use. Suitable for stainless steel, ceramics, glass, plastic and soft metal surfaces. The pH neutral, water based formulation is easily rinsed off in manual pre-wash applications.



This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at getinge.com/contact

www.getinge.com

GETINGE