

Getinge Clean Rinse Aid

Speed up drying times

Bio-compatible additive used during the final rinse to help drying

Getinge Clean Rinse Aid is a surfactant based rinse aid for use in washer-disinfectors to ease drying and achieve spot free rinsing of medical devices.

Suitable for use on stainless steel, ceramics, glass, plastics and resistant soft metal surfaces. Rinse Aid will not interfere with steam or ethylene oxide sterilization processes.



Available in 5L, 10L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at getinge.com/contact

www.getinge.com

GETINGE 