

# Getinge Clean Universal Detergent

Mild alkaline based detergent

Getinge Clean Universal Detergent is suitable for most types of materials

Getinge Clean Universal Detergent is a mild alkaline detergent for use in washer-disinfectors, ultrasonic baths and manual cleaning of medical devices.

It is a low foaming, non-abrasive, instrument cleaning detergent that rapidly removes, dissolves and disperses most types of soil.

Suitable for use on alkali resistant surfaces such as stainless steel, ceramics, glass, resistant plastics and soft metals.



Available in 5L, 10L, 25L and 205L containers

This document is intended to provide information to an international audience outside of the US.

For more detailed information, please contact your local Getinge representative at [getinge.com/contact](http://getinge.com/contact)

[www.getinge.com](http://www.getinge.com)

**GETINGE** 