

# Getinge Pack Standard Packaging

## Sterilization rolls and pouches

Providing high quality and strong microbial barriers for all your processed hospital instruments.

Sterilization rolls and pouches provide healthcare workers an easy and safe method to enclose medical devices intended to be sterilized in steam, ethylene oxide (EO) or formaldehyde sterilizers.

The rolls and pouches, made of high grade medical paper and laminated multilayer PET/PP copolymer film, ensure an effective microbial barrier for all your processed hospital instruments. Water based and non-toxic process indicators printed on the medical paper surface provide an accurate visual confirmation by color change and will help to differentiate between processed and unprocessed packages.

A wide range of standard sizes allows you for optimum choice of a correct-sized pack for each instrument and for double packing purposes.

Getinge Pack sterilization rolls and pouches are CE marked and in compliance with international standards ISO 11607-1, ISO 11140-1 and EN 868-5.

This document is intended to provide information to an international audience outside of the US.

