



Linear IAB Catheter

7.5 Fr. IAB Technical Specifications

Technical Specifications

Linear 7.5 Fr. IAB Catheter with Insertion Kit

Part Number	Balloon					
	Nominal Volume (cc)	Nominal Length (mm)	Balloon Dia. at 2 psi (mm)	Insertable Length (mm)	Pressure Lumen ID Size mm (in.)	Guidewire Size mm (in.)
0684-00-0480-01	40	258	15.0	723	0.68 mm (0.027")	0.64 mm (0.025")
0684-00-0479-01	34	221	15.0	723	0.68 mm (0.027")	0.64 mm (0.025")
0684-00-0478-01	25	165	15.0	723	0.68 mm (0.027")	0.64 mm (0.025")



Ordering Information



Product	Part Number
Linear 7.5 Fr. 40 cc IAB with Insertion Kit	0684-00-0480-01
Linear 7.5 Fr 40 cc IAB with Insertion Kit and Arrow Pump Adapter (APA)	0684-00-0480-02
Linear 7.5 Fr. 34 cc IAB with Insertion Kit	0684-00-0479-01
Linear 7.5 Fr. 34 cc IAB with Insertion Kit and Arrow Pump Adapter (APA)	0684-00-0479-02
Linear 7.5 Fr. 25 cc IAB with Insertion Kit	0684-00-0478-01
Linear 7.5 Fr 25cc IAB with Insertion Kit and Arrow Pump Adapter (APA)	0684-00-0478-02
Linear 7.5 Fr. 40 cc and 34 cc Insertion Kit includes:	0884-00-0019-13
<ul style="list-style-type: none"> • One 15 cm (6") Reinforced Introducer Sheath with hemostasis valve • One Introducer Dilator • One 1.3 mm (18G, 0.050") Angiographic Needle • One Vessel Dilator • One 0.64 mm (0.025") X 145 cm, 3 mm J, PTFE coated, Stainless Steel Guidewire • One 183 cm (6') Catheter Extender Tubing (40 cc and 34 cc only) • One Luer Cap • One Three-Way Stopcock • Two 122 cm (4') Pressure Tubings 	
Linear 7.5 Fr. 25 cc Insertion Kit includes:	0884-00-0019-12
All above contents except 6' Catheter Extender Tubing replaced with 5' Catheter Extender Tubing	
Arrow Pump Adapter (APA) 5-Pack 25 cc to 40 cc	0684-00-0510-01
Statlock for IAB (20 pack) – Catheter Staibilization Device	0684-00-0472

Configuration

Packaging	One sterile pouch, corner peel
Cycle Life	One procedure
IAB Shelf Life	3 years from date of manufacturing

For more product information, contact your local representative.

MCA00000693 REVC · Getinge, and **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries · StatLock is a registered trademark of C.R. Bard, Inc. or an affiliate. Copyright 2024 Datascope Corp. · All rights reserved ·  Refer to Instructions for Use for current indications, warnings, contraindications and precautions · 10/2024

 Getinge · 45 Barbour Pond Drive · Wayne, NJ 07470 · USA
 Datascope Corp. · 15 Law Drive · Fairfield, NJ 07004 · USA

www.getinge.com

GETINGE 