



Maquet Moduevo Ceiling Supply Units

Accessories catalog

This document is intended to provide information to an international audience outside of the US.

GETINGE 



Maquet Moduevo accessories

– A wealth of options

With Maquet Moduevo Ceiling Supply Units, clinicians can spend less time accessing or moving equipment, so they have more time for patient care.

Maquet Moduevo has been designed with a singular purpose in mind: to enhance provider-patient interactions with safety, efficiency and long-term ownership considerations.

Maquet Moduevo does so with extraordinary simplicity, delivering effortless maneuverability in the Operating Room, Intensive Care Units, Emergency Departments, and other space-constrained clinical environments – with unprecedented ease in customization, and an uncomplicated approach to future-proof flexibility.

In addition, Maquet Moduevo Ceiling Supply Units save space, and improve safety by providing effective cable management.

They also ensure vital utilities and equipment are readily accessible, streamlining workflows and raising efficiency.

A broad choice of accessories further enhances their versatility.





About this catalog

This catalog provides an overview of the Maquet Moduevo range of ceiling supply units. It is divided into **four categories** based on the products' function: positioning, holding, illuminating and storing.

Position accessories include a range of handles with various brakes and positioning options. **Hold** items comprise shelves, drawers, monitor arms, rails, adaptors and infusion equipment (e.g. IV poles).

The **illuminate** category presents lighting options for diverse situations and needs. And **storage** accessories encompass cable management products such as boxes, plus holders and baskets.

Each of these four main categories is subdivided into sections dedicated to specific types of items, such as those listed above. These sections provide readers with more detailed information and technical specifications for individual accessories.

Maquet Moduevo accessories

– Content

To position	6
Handles	7
To hold	9
Shelves	10
Cabinets/drawer units	13
Flat screen monitor arms	15
Rails	18
Infusion equipment	20
Adaptors/docking systems	24
To illuminate	27
Ambient lighting	27
To store	29
Cable management solutions	30
Baskets/catheter holders	32
Easy clamp holders	35

To position

– Maquet Moduevo accessories

A key feature of Maquet Moduevo is their outstanding maneuverability. Staff can quickly and precisely position the ceiling supply units where they are needed.

To this end, there are a variety of handles and control devices available. For example, handles may be attached horizontally or vertically, and to various sides of the ceiling supply unit. Positioning can be done simply mechanically (by hand), or it can involve more advanced technology –

e.g. ceiling supply unit height can be adjusted via the push of a button, and brakes released just by grasping a capacitive sensing handle. Whatever the scenario, Getinge offers ergonomic solutions that are easy for staff to use, helping them respond rapidly to changing situations.



Handles

Handles serve a simple yet vital purpose: they enable staff to position the ceiling supply unit easily and precisely. To this end, a variety of designs are available for Maquet Moduevo – including vertical and horizontal handles, located on the front or back of the unit in line with customers’ specific needs.

Maquet Moduevo handles can incorporate innovative Grab & Move Capacitive Sensing Technology: the user simply has to grasp the handle with one hand to release the brakes. The ceiling supply unit can then be moved to the desired position with minimal physical effort. In addition, handles are compatible with various Maquet Moduevo arms or beams – including the motorized, height-adjustable Energy and Cargolift.



Horizontal control handle



Vertical control handle



Features

- Can be directly mounted on EPoS System
- Includes Grab & Move Capacitive Sensing Technology if brakes option has been selected
- Enables flexible positioning
- Rounded edges are easy to clean

Technical information: horizontal control handle	Order number
With 1-axis or 2-axis brake control	SU000600054
Technical information: vertical control handle	Order number
Without brake control	SU000600059
With 1-axis or 2-axis brake control	SU000600061

Horizontal control handle (for height adjustment)



Features

- Can be directly mounted on EPoS System
- Includes Grab & Move Capacitive Sensing Technology if brakes option has been selected to release brakes and maneuver ceiling supply unit
- Enables flexible positioning
- Rounded edges are easy to clean
- For use with Energy and Cargolift
- Ceiling supply unit height can be adjusted at the push of a button

Technical information: horizontal control handle	Order number
Without brake control	SU000600055
With 1-axis or 2-axis brake control	SU000600057
Optional 200 mm / 7.87" spacer for horizontal handle	SU000600058

Vertical control handle (for height adjustment)



Technical information: vertical control handle	Order number
Without brake control	SU000600062
With 1-axis or 2-axis brake control	SU000600064

Shelf handle



Features

- Directly mounted on shelves to save space
- Various options for brake control

Technical information	Order number
Without brake control for height adjustment W = 485 mm / 19.0" W = 685 mm / 27.0"	SU000600067 SU000300042
With 1-axis or 2-axis brake control for height adjustment W = 485 mm / 19.0" W = 685 mm / 27.0"	SU000600069 SU000300044
With 1-axis or 2-axis brake control for standard rotation W = 485 mm / 19.0" W = 685 mm / 27.0"	SU000600066 SU000300041

Ceiling supply unit control device



Features

- Wired control device with an electrical clamp
- Maneuvers the ceiling supply unit up and down
- Clamp can be mounted on EPoS System

Technical information	Order number
For Energy and Cargolift Ceiling Supply Units	SU000600070

To hold

– Maquet Moduevo accessories

A key attraction of ceiling supply units, including Maquet Moduevo, is that a variety of equipment can be added to a central, easily moved station.

The low voltage power supplied by the EPoS System paired with a broad portfolio of accessories mean hospital staff can quickly and flexibly configure the ideal environment for their task at hand.


With this in mind, the to hold category comprises the largest range of accessories for Maquet Moduevo: shelves, drawers and cabinet units, flat screen monitor arms, rails, infusion equipment, plus adapters and docking systems. They all play a key role in ensuring critical equipment is readily accessible whenever it is required.



Shelves

Shelves are seemingly simple, yet fulfill an important purpose by ensuring critical equipment is on hand for direct access and use. By organizing items and cutting out clutter, they make better use of space, and enhance efficiency.

Moreover, they feature high load capacity and are exceptionally durable. Shelves for Maquet Moduevo Units can be attached to EPoS System via mechanical or electric clamps – the latter maintains a low-voltage power connection.



Shelf with mechanical clamp and two rails



Features

- Available in three sizes, compatible with the shelf extension
- Load capacity: 30 kg / 66 lbs with single clamp, 50 kg / 110.2 lbs with two clamps
- Two side rails (25 x 10 mm / 0.98" x 0.39"); load capacity: 10 kg / 22 lbs each
- Optional mechanical handle
- All single clamp shelves can be upgraded to double clamp shelves and vice versa

Technical information	Order number
With single clamp and no handle for distributors W x D = 485 mm x 450 mm / 19.0" x 17.7" W x D = 685 mm x 450 mm / 27.0" x 17.7"	SU000600083 SU000600085
With single clamp and mechanical handle for distributors W x D = 485 mm x 450 mm / 19.0" x 17.7" W x D = 685 mm x 450 mm / 27.0" x 17.7"	SU000600088 SU000600092
With double clamps and no handle for distributors W x D = 485 mm x 450 mm / 19.0" x 17.7" W x D = 685 mm x 450 mm / 27.0" x 17.7" W x D = 885 mm x 450 mm / 34.9" x 17.7"	SU000600084 SU000600086 SU000600087
With double clamps and mechanical handle for distributors W x D = 485 mm x 450 mm / 19.0" x 17.7" W x D = 685 mm x 450 mm / 27.0" x 17.7"	SU000600089 SU000600093
With double clamps and no handle for racks W x D = 485 mm x 450 mm / 19.0" x 17.7" W x D = 685 mm x 450 mm / 27.0" x 17.7" W x D = 885 mm x 450 mm / 34.9" x 17.7"	SU000300066 SU000300069 SU000300072

With double clamps and mechanical handle for racks

W x D = 485 mm x 450 mm / 19.0" x 17.7"
 W x D = 685 mm x 450 mm / 27.0" x 17.7"
 W x D = 885 mm x 450 mm / 34.9" x 17.7"

SU000300067
 SU000300070
 SU000300073

Single mechanical shelf clamp for upgrade (right side)

SU000600104

Double clamp shelf spacer 75 mm / 2.95",
 50 kg / 110.2 lbs loading

SU000600105

Single clamp shelf spacer 75 mm / 2.95", 30 kg / 66 lbs loading

SU000600130

Shelf with handle and single electrical clamp



Features

- Available in two sizes
- Compatible with the shelf extension
- Load capacity (with single clamp): 30 kg / 66 lbs
- Load capacity (with two clamps): 50 kg / 110.2 lbs
- Two side rails (25 mm x 10 mm / 0.98" x 0.39"); load capacity: 10 kg / 22 lbs each
- Optional Grab & Move Handle for easy one-handed positioning

Technical information

Order number

With single electrical clamp and two rails for distributors

W x D = 485 mm x 450 mm / 19.0" x 17.7"
 W x D = 685 mm x 450 mm / 27.0" x 17.7"

SU000600090
 SU000600094

With one electrical clamp, one mechanical clamp and two rails for distributors

W x D = 485 mm x 450 mm / 19.0" x 17.7"
 W x D = 685 mm x 450 mm / 27.0" x 17.7"

SU000600091
 SU000600095

With one electrical clamp, one mechanical clamp and two rails for racks

W x D = 485 mm x 450 mm / 19.0" x 17.7"
 W x D = 685 mm x 450 mm / 27.0" x 17.7"
 W x D = 885 mm x 450 mm / 34.9" x 17.7"

SU000300068
 SU000300071
 SU000300074

Sliding keyboard shelf



Features

- Smart storage for keyboards; folds away when not in use
- Load capacity: 5 kg / 11 lbs

Technical information

Order number

W x D = 485 mm x 450 mm / 19.0" x 17.7"

SU000600100

Shelf extension



Features

- No need for special tools, simple installation of extension plates onto side rails
- Provides more space on the shelf to accommodate larger or additional equipment
- Compatible with all shelves

Technical information

Order number

100 mm / 3.94" wide for each extension plate

SU000600101

Anti-sliding bars



Features

- Available in three sizes
- Secures the positioning of devices and accessories on Moduevo shelves

Technical information	Order number
For shelf 485 mm / 19.0" wide	SU000600151
For shelf 685 mm / 27.0" wide	SU000600152
For shelf 885 mm / 34.9" wide	SU000600153

Swivel shelf



Features


- Rotary arm allows shelf to be extended outwards
- Range of movement can be further increased by means of the optional double-articulation arm
- Smart storage: shelf can be folded away when not in use

Technical information	Order number
W x D = 400 mm x 360 mm / 15.7" x 14.2", load capacity: 15 kg / 33 lbs	
Single arm (L = 200 mm / 7.9")	SU000600033
Double arm (L = 300 mm + 300 mm / 11.8" + 11.8")	SU000600034

Cabinets and drawer units

Drawers can hold a variety of essential items, making them useful across a range of medical disciplines and healthcare environments. Various cabinet units with one or two drawers can be easily mounted onto the ceiling supply unit's EPoS System.

They come with an in-drawer light as standard, helping hospital staff locate instruments and materials safely and efficiently in dark environments.



Cabinet unit with one drawer



Features

- Available in two sizes
- Can be mounted below any type of shelf
- Total load capacity: 10 kg / 22 lbs
- Depth: 100 mm / 3.9"
- In-drawer ambient light
- Optional drawer lock (SU000600102)
- Detachable and interchangeable, facilitating the transfer of contents between stations and rooms

Technical information	Order number
W x D x H = 485 mm x 330 mm x 100 mm / 19.0" x 13.0" x 3.9"	SU000600096
W x D x H = 685 mm x 330 mm x 100 mm / 27.0" x 13.0" x 3.9"	SU000600097

Cabinet unit with two drawers



Technical information	Order number
W x D x H = 485 mm x 330 mm x 100 mm / 19.0" x 13.0" x 3.9" (per unit)	SU000600098
W x D x H = 685 mm x 330 mm x 100 mm / 27.0" x 13.0" x 3.9" (per unit)	SU000600099

Cabinet unit with one high drawer



Features

- Available in two sizes
- Can be mounted below any type of shelf
- Total loading capacity: 10 kg / 22 lbs
- Depth: 250 mm / 9.8"
- In-drawer ambient light
- Optional drawer lock (SU000600102)
- Detachable and interchangeable, facilitating the transfer of contents between stations and rooms

Technical information	Order number
W x D x H = 485 mm x 330 mm x 250 mm / 19.0" x 13.0" x 9.8"	SU000600149
W x D x H = 685 mm x 330 mm x 250 mm / 27.0" x 13.0" x 9.8"	SU000600150

Cabinet unit with drawer and e-lock



Features

- Prevent unauthorized access to objects stored in the drawer
- RFID technology
- Powered by EPoS with 24 V DC
- The authorization of user cards can be programmed
- Compatible with customers' RFID cards with transmission frequency of 13.56 MHz
- Compatible with Moduevo shelves (width of 485 mm / 19.0") with electrical clamp
- Depth: 100 mm / 3.9"
- In-drawer ambient light

Technical information	Order number
Single drawer with e-lock – 485 mm / 19.0"	SU000600155
Double drawer with one e-lock – 485 mm / 19.0"	SU000600156
Double drawer with two e-locks – 485 mm / 19.0"	SU000600157
Electrical lock card	SU000600154

Flat screen monitor arms

In ORs, ICUs and beyond, flat screen monitors provide critical information and visuals at a glance. Maquet Moduevo's range of monitor arms, or holders, enables displays to be moved and tilted – so they are always clearly visible, and exactly where they are required.

Moreover, the arms feature a sleek design, and are easy to adjust – for smaller screens up to larger, 14 kg/ 30.9 lbs monitors. Their wide range of movement contributes to an ergonomic working environment.



Flat screen monitor holder



Features

- Holds flat screen with VESA 75/100 adaptor

Technical information	Order number
Simple holder for EPoS	SU001900050

Tilt and rotation monitor holder



Features

- Holds flat screen with VESA 75/100 adaptor
- Load capacity: 14 kg/ 30.9 lbs

Technical information	Order number
L = 100 mm/3.94"	SU000600018

Fixed-height single arm monitor holder



Features

- Available in two lengths
- Includes swivel arm and tiltable screen holder
- 180° side-to-side rotation
- Load capacity: 14 kg / 30.9 lbs
- VESA 75/100 adaptor included
- Cables are routed through the arm to guarantee a safe and hygienic workplace

Technical information	Order number
L = 300 mm / 11.8"	SU000600019
L = 400 mm / 15.7"	SU000600020

Fixed-height double arm monitor holder



Features

- Available in three lengths
- Includes swivel arm and tiltable screen holder
- 180° side-to-side rotation
- Load capacity: 14 kg / 30.9 lbs
- VESA 75/100 adaptor included
- Cables are routed through the arm to guarantee a safe and hygienic workplace

Technical information	Order number
L = 300 mm + 300 mm / 11.8" + 11.8"	SU000600021
L = 400 mm + 300 mm / 15.7" + 11.8"	SU000600022
L = 400 mm + 400 mm / 15.7" + 15.7"	SU000600023

Height-adjustable monitor holder



Features

- Choice of single- or double-articulation arm
- Depth, angle and height adjustable
- Side-to-side rotation
- Vertical adjustment
- Compatible with monitors up to 14 kg / 30.9 lbs
- Range of movement: 700 mm / 27.5"
- VESA 75/100 adaptor

Technical information	Order number
L = 300 mm / 11.8"	SU000600024
L = 400 mm + 300 mm / 15.7" + 11.8"	SU000600025

Keyboard holder



Features

- Easy to install and use
- Equipped with mouse pad
- Load capacity: 5 kg / 11 lbs

Technical information

Keyboard holder

Order number

SU000600026

Down post for monitor arms



Features

- For monitors mounted on non-height-adjustable and/or height-adjustable monitor holders

Technical information

For non-height-adjustable monitor holder (L = 196 mm / 7.7")

For height-adjustable monitor holder (L = 196 mm / 7.7")

Universal adapter (load capacity: 14 kg / 30.9 lbs)

Order number

SU000600027

SU000600028

SU000600029

Rails

Rail systems for Maquet Moduevo Ceiling Supply Units can be easily installed on the EPoS System – there is no need for specialized tools.

These rails are available in a number of lengths and clamps, meaning they can be employed for various Maquet Moduevo Distributors. They allow diverse devices and accessories, such as baskets, flow meters and glove boxes to be mounted on the units, increasing their versatility and saving space.



Standard rail



Features

- Available in four sizes for various panels and Maquet Moduevo Distributors
- Simple installation on EPoS System
- Load capacity: 30 kg / 66 lbs

Technical information	Order number
200 mm / 7.9" with single mechanical clamp	SU000600008
400 mm / 15.8" with single mechanical clamp	SU000600009
400 mm / 15.8" with double mechanical clamp for back of Mezzo Distributors	SU000600010
400 mm / 15.8" with double mechanical clamp for front of Slim / Mezzo / Twin Distributors	SU000600011

Wrap-around rail



Features

- Improves efficiency by enabling staff to attach equipment anywhere around the distributor
- Available in three sizes for various types of distributors
- Load capacity: 15 kg/33 lbs

Technical information	Order number
For Slim Distributor	SU000600129
For Mezzo Distributor	SU000600016
For Twin Distributor	SU000600015

Half-wrap rail



Features

- Improves efficiency by enabling staff to attach equipment across half of the distributor
- Available in two sizes for various types of distributors
- Load capacity: 15 kg/33 lbs

Technical information	Order number
For Twin Distributor	SU000600013
For Mezzo Distributor	SU000600014

Swivel rail



Features

- Extendible for broader range of movement
- Can be folded away when not in use
- Simple installation without specialized tools
- Load capacity: 3 kg/6.6 lbs

Technical information	Order number
Swivel rail (max. L = 390 mm/15.4")	SU000600012

Infusion equipment

Infusions are part and parcel of healthcare – so making them more efficient can have a significant overall impact.

Maquet Moduevo IV Poles can accommodate up to 50 kg / 110 lbs of infusion pumps, bags and more. Plus, the pole height and position of the arm attached to the ceiling supply unit can be adjusted and securely locked in place. And the Maquet M-Shift IV Transfer System enables IV equipment to be moved quickly yet safely, easing nurses' workloads.



Standard IV pole



Features

- Flexible solution – various ranges of movement with single- or double-articulation arms; arms are rotary
- Each arm can be locked into position thanks to the smart and ergonomic “one-finger knob”
- Simple installation on EPoS System

Technical information: standard IV pole	Order number
One standard IV pole with four hooks in a straight line, Ø 25 mm / 1.0", single clamp, load capacity: 20 kg / 44 lbs	
Single-articulation arm (L = 180 mm / 7.1")	SU000600001
Double-articulation arm (L = 180 mm + 180 mm / 7.1" + 7.1")	SU000600002
Single-articulation arm (L = 300 mm / 11.8")	SU000600141
Double-articulation arms (L = 300 mm + 300 mm / 11.8" + 11.8")	SU000600142
With optional IV pole	SU000600007

Height-adjustable IV pole



Features

- Flexible solution – various ranges of movement with single- or double-articulation arms; arms are rotary
- Each arm can be locked into position thanks to the smart and ergonomic "one-finger knob"
- Simple installation on EPoS System
- Pole height can be easily adjusted and locked in place

Technical information: height-adjustable IV pole	Order number
One height-adjustable IV pole with four hooks in a cluster, Ø 25 mm / 1.0", single clamp, load capacity: 20 kg / 44 lbs	
Single-articulation arm (L = 180 mm / 7.1")	SU000600003
Double-articulation arm (L = 180 mm + 180 mm / 7.1" + 7.1")	SU000600004
Single-articulation arm (L = 300 mm / 11.8")	SU000600143
Double-articulation arm (L = 300 mm + 300 mm / 11.8" + 11.8")	SU000600144
With optional IV pole	SU000600007

Strong IV pole



Features

- Pair of arms mounted on EPoS
- Flexible solution – various ranges of movement with single- or double-articulation arms; arms are rotary
- High load capacity

Technical information	Order number
One strong IV pole with four hooks in a straight line, Ø 25 mm / 1.0", double clamps, load capacity: 50 kg / 110 lbs	
Pair of single-articulation arms (L = 180 mm / 7.1"), 50 kg / 110 lbs	SU000600005
Pair of double-articulation arms (L = 180 mm + 180 mm / 7.1" + 7.1"), 50 kg / 110 lbs	SU000600006
Single-articulation arm (L = 300 mm / 11.8"), 40 kg / 88 lbs	SU000600145
Pair of double-articulation arms (L = 300 mm + 300 mm / 11.8" + 11.8"), 50 kg / 110 lbs	SU000600146

Maquet M-Shift



Features

- Flexible infusion transfer system: two angle- and height-adjustable IV poles
- Net load capacity of 50 kg/110 lbs for pumps and bags
- Small footprint – easy to fit into elevators
- Only one nurse is needed to safely transfer the infusion system from ceiling supply unit to bed and vice versa

Technical information	Order number
Maquet M-Shift	SU0005990010
Maquet M-Shift Bed Adaptor Set connects M-Shift with ArjoHuntleigh ICU Bed; includes bed pole and adaptor E	SU0005990013
Adaptor E for ArjoHuntleigh ICU Bed to hold docking pole	SU0005060002
Bed pole, can be combined with Maquet M-Shift; includes stopper; requires bed adaptor	SU0005060005
Maquet M-Shift Trolley to carry Maquet M-Shift Clamp System; includes trolley, one standard rail, one basket, one bed pole, one cable management	SU0005060006

Maquet M-Shift Docking System



Features

- Enables Maquet M-Shift Docking Pole to be fixed to Maquet Moduevo Ceiling Supply Units

Technical information	Order number
Adaptor for ArjoHuntleigh, Hillroom and Stryker Beds, subject to bed compatibility validation. Includes bed adaptor and docking pole	
Docking pole for Maquet Moduevo 800 Modules	SU000600032



88% 88:88 88

PHILIPS CODE: 4830

HAQUET

GERINGE



Adaptors / docking equipment

GCX Adaptor



Features

- Accommodates GCX® accessories on the EPoS System

Technical information	Order number
Adaptor for GCX Monitor Holders (loading capacity: 14 kg / 30 lbs)	SU000600031
Elongated GCX Adaptor on EPoS (loading capacity: 20 kg / 44 lbs)	SU000600128
Adaptor on suspension tube for monitor arm and GCX Adaptor	SU000600030

Beam Plus Adaptor



Features

- Integrated adaptor installed below the primary ceiling supply unit axis
- Used to mount examination lights, surgical lights, and monitor or camera holders
- Adaptor and ceiling supply unit beam rotate independently

Technical information	Order number
For surgical light Maquet Lucea 40 / 50 / 100 and Maquet Volista Access 40 SF / DF For anchor I For anchor II	SU000300046 SU000900010
For OR lights Maquet Lucea 50 / 100, SAT XS 0912 / 0913 / 1212 / 1213, SAT XO 12 / 13, SAT Orchide 12 / 13, EQTFHS010 SAT12 / 13, EQTMHS013 SAT 12 / 13, EQTMHS021 SAT12 / 13, EQTMHD231 SAT 12 / 13, EQTXHS010 SAT12 / 13, EQTXHS021 SAT12/13 For anchor I For anchor II	SU000300047 SU000900011
BPA reservation	SU000900014

Bridge Plus Adaptor



Features

- Integrated adaptor installed above Bridge classic or advance beam
- Used to mount OR and examination lights

Technical information	Order number
For OR lights Maquet Lucea 50 /100 Bridge Plus SAT	SU000200078
For OR light Maquet Lucea 40, 50 SF / DF, 100 SF / DF, Maquet Volista Access 40 SF / DF and 60 SF / DF Bridge Plus SB	SU000200079

Anesthesia docking system



Features

- Docking system for anesthesia machines; designed to streamline workflows and improve safety
- High load capacity; can accommodate diverse anesthesia machines

Technical information	Order number
For anesthesia machines and MIS Trolley on Energy, only compatible with Twin 406 / 606 / 806, loading capacity: 180 kg / 396.8 lbs	SU000300026

Ventilator docking holder



Features

- Docking holder for Getinge ventilators
- Designed to streamline workflows and improve safety

Technical information	Order number
For mounting Servo-u/-n ventilators with bed / pendant bracket on EPoS, loading capacity: 30 kg / 66 lbs	SU000600132

Gas cylinder holder



Features

- Docking holder for gas cylinders
- Flexible design of support plate at bottom fits cylinders with a diameter from 100 mm to 180 mm / 3.9" to 7.1"
- Basic holder for gas cylinders no higher than 450 mm / 17.7"
- Support clamp adds stability for gas cylinders with height from 450 mm to 865 mm / 17.7" to 34.1"
- Loading capacity: 30 kg / 66 lbs

Technical information	Order number
Basic gas cylinder holder	SU000600147
Support clamp for high gas cylinder	SU000600148

OR tables charging base solutions



Features

- Charging base for Getinge Operating Table Hand Control (1009.25A0) that can be continuously powered by EPoS with 24 V DC
- Compact design
- Compatible with Maquet Otesus and Maquet Meera Operating Tables

Technical information	Order number
Operating table remote control charging base on EPoS	SU000600131

To illuminate

– Maquet Moduevo accessories

Better visibility is key to better working environments and streamlined tasks. A soft glow contributes to a soothing atmosphere for convalescing patients, and supports natural circadian cycles.

Meanwhile, staff can locate items, complete tasks and continue treatment without disruption. Against this backdrop, a system of ambient lights for Maquet Moduevo enhances efficiency by effectively illuminating dark ORs and ICUs.

This category comprises a range of ambient lights that can be used with various beams and Maquet Moduevo Distributors – and be placed on the top, bottom or side of the ceiling supply unit to eliminate shadows. Additionally, brightness can be adjusted with a single control unit, and dimmed down to just 16%.



Ambient light

Special ambient lighting does away with harsh glare and shadows. Instead, the soft illumination creates a soothing environment ideal for resting patients – without compromising visibility for staff as they go about their tasks.

The lights are positioned on top of the ceiling supply unit beam, on the distributor's sides, and at the bottom of the ceiling supply unit, for optimal results and improved safety.

Features

- For all types of ambient lights; just one dimmer control (SU000400001) is required, and can be used to control all lights
- Certain ceiling supply units might not allow all the light placement options – please consult your Getinge representative for advice
- Brightness can be adjusted from 100% to 16% to meet specific needs

Technical information	Order number
Top illumination for upper beam	SU000600076
Side illumination for Slim and Mezzo Distributors	SU000600077
Side illumination for Twin Distributors	SU000600078
Bottom illumination for Slim, Mezzo and Twin Distributors	SU000600081
Bottom illumination for Sky Distributors	SU000300053
Bottom illumination for Moduevo Screen	SU001900013

Moduevo Bridge ambient light



Features

- Various regular lighting solutions on Moduevo Bridge, offering soft illumination for better working atmosphere
- Reading / task light provides two different modes with illumination level less than 500 lux for the reading and 5,000 lux for the task light
- Somnus Light provides innovative daylight simulation and with the same optical system it provides also ambient light with six different colors

Technical information	Order number
Night light	SU000400085
Reading / task light	SU000400084
Side illumination for Slim and Mezzo Distributors on Bridge	SU000400086
Bottom illumination for Slim and Mezzo Distributors on Bridge	SU000400087
Bottom illumination for Sky Distributors on Bridge	SU000300086
Somnus and colored ambient lights on Bridge horizontal beam	SU000400088

Ambient light keypads



Features

- Touch panel
- Smooth illumination adjustment
- Hygienic design for easy cleaning
- Flexible solution enables users to control ambient lights from either the pendant or the wall

Technical information	Order number
Ambient light dimmer controls for regular Moduevo CSU	
Dimmer control on distributors	SU000400001
Dimmer control on Sky classic	SU000300076
A pair of dimmer controls on distributors and on the wall	SU001900014
Ambient light dimmer controls for Moduevo Bridge	
Basic dimmer control for Bridge	SU000400089
Premium dimmer control for Bridge	SU000400090
Wall mounted basic dimmer control for Bridge	SU000400091
Wall mounted premium dimmer control for Bridge	SU000400092
A pair of basic dimmer controls on Bridge and on the wall	SU000400094
A pair of premium dimmer controls on Bridge and on the wall	SU000400095

To store

– Maquet Moduevo accessories

The storage category of Maquet Moduevo accessories is dedicated to simple, yet effective solutions for managing cables, and storing documents, gloves, equipment such as catheters, and more.

In line with the general philosophy behind ceiling supply units, items are kept off the ground. This improves access to the patients – and creates a safer and more hygienic workplace for staff.

And this degree of organization has additional benefits. Because every item has its place, staff know exactly where to find what they are looking for. This makes for more efficient workflows, and allows healthcare professionals to turn their attention to more important matters.



Cable management solutions

Conventional ORs and ICUs might have several gas tubes and electric cables snaking across the floor. Not only do these exposed wires add to the surface area that has to be disinfected – they are also tripping hazards.

This is where Maquet Moduevo Ceiling Supply Units come into play: by incorporating smart cable management solutions, they keep tubes off the floor – for cleaner, safer healthcare environments.



Cable winding wheel



Features

- Cable management provides protection from dust, simplifying cleaning and contributing to a hygienic environment
- Cables can be stored in a tidy, organized way
- Four options available; to be mounted on EPOS, below shelves, and using a small hook

Technical information	Order number
Cable winding wheel with clamp	SU000600035
Cable winding wheel for below shelves	SU000600103

Cable holders



Features

- Ideal solution for gas tubes and electrical cords
- Improves hygiene and safety

Technical information	Order number
Rack cable holder, set of 5 units	SU000300065
Cable holder for rails, set of 10 units	SU000600017

Cable box



Features

- Ideal solution for gas tubes and electrical cords
- Improves hygiene and safety

Technical information	Order number
For installation on open channel of Twin Distributor	SU000600036

Baskets /catheter holders

Storage baskets come in a variety of shapes and sizes. There are, for example, multiple options for holding catheters, important documents, gloves and more – ensuring these items are organized and easy to find. The baskets can be hooked or installed directly onto rails, and are made from hygienic stainless steel.



Catheter holder



Features

- Simple installation: two hooks for attachment onto rails
- Hygienic: stainless steel

Technical information

H = 500 mm / 19.7"

Order number

ARD513768999

Small basket



Features

- Simple installation: two hooks for attachment onto rails
- Hygiene: stainless steel


Technical information

W x D x H = 170 mm x 130 mm x 85 mm / 6.7" x 5.12" x 3.35"

Order number

SU000600126

Catheter holder



Features

- Basket for catheter, available in two sizes
- Hygiene: stainless steel

Technical information	Order number
Small: W x D x H = 90 mm x 90 mm x 350 mm / 3.54" x 3.54" x 13.8"	SU000600125
Large: W x D x H = 150 mm x 100 mm x 480 mm / 5.9" x 3.9" x 18.9"	SU0005990000

Wire cylinder holder



Features

- Jar holder
- Simple installation: direct attachment onto rail
- Hygiene: stainless steel

Technical information	Order number
Ø 110 mm / 4.3"	ARD513769999

Wire document holder



Features

- Simple installation: direct attachment onto rail
- Hygiene: stainless steel

Technical information	Order number
Wire document holder	ARD686400051

Bag holder



Features

- Simple installation
- Hygiene: stainless steel

Technical information	Order number
W x D = 250 mm x 170 mm / 9.9" x 6.7"	ARD517052914

Gloves holder



Features

- Simple installation: direct attachment onto rail
- Hygiene: stainless steel

Technical information	Order number
H = 200 mm / 7.9"	ARD513766999

Easy clamp holders

In clinical environment, small equipment needs to be placed close to patients to streamline clinical workflows.

Easy clamp holders, based on snap action clamps, can be quickly attached to or detached from standard DIN rail without any tool, thus can easily meet unique needs for well organization of equipment and bedside space management.



Glove box holder



Features

- Simple installation on DIN rail with easy clamp
- Hygiene: stainless steel

Technical information

W x D x H = 135 mm x 110 mm x 264 mm / 5.3" x 4.3" x 10.4"

Order number

SU000600136

Ambu-bag / resuscitator holder



Features

- Simple installation on DIN rail with easy clamp

Technical information

For mounting ambu-bag / resuscitator

Order number

SU000600134

Flexible wipes canister holder



Features

- Simple installation on DIN rail with easy clamp
- Flexible design for wipes canisters with different diameters

Technical information	Order number
For mounting wipes canisters with diameter from 98 mm to 150 mm / 3.8" to 5.9"	SU000600133

Scanner holder



Features

- Simple installation on DIN rail with easy clamp

Technical information	Order number
Inner dimension W x D x H = 80 x 102 x 50 mm / 3.1" x 4.0" x 2.0"	SU000600138

Welch Allyn Aneroid Holder



Features

- Simple installation on DIN rail with easy clamp

Technical information	Order number
For mounting Welch Allyn 767 Series Aneroid	SU000600137

Welch Allyn Ophthalmoscope Holder



Features

- Simple installation on DIN rail with easy clamp

Technical information	Order number
For mounting Welch Allyn Green Series 777 Ophthalmoscope	SU000600139

Welch Allyn Ophthalmoscope Tip Container Holder



Features

- Simple installation on DIN rail with easy clamp

Technical information	Order number
For mounting Welch Allyn Ophthalmoscope Tip Container	SU000600135

Welch Allyn Thermometer Holder



Features

- Simple installation on DIN rail with easy clamp

Technical information	Order number
For mounting Welch Allyn Suretemp 690/692 Thermometer	SU000600140



With a firm belief that every person and community should have access to the best possible care, Getinge provides hospitals and life science institutions with products and solutions aiming to improve clinical results and optimize workflows. The offering includes products and solutions for intensive care, cardiovascular procedures, operating rooms, sterile reprocessing and life science. Getinge employs over 10,000 people worldwide and the products are sold in more than 135 countries.

This document is intended to provide information to an international audience outside of the US.

Maquet (Suzhou) Co., Ltd. · No.158 Fangzhou Road · Suzhou Industrial Park · 215024 Suzhou · China · +86 (0) 512 62839880

www.getinge.com