

GET250, GET350, GET750

Blanket Warming Cabinets

Single chamber blanket warming cabinets. Interior surface temperatures are monitored and balance heat throughout the cabinets. The chambers have an adjustable temperature range (32 °C to 71 °C/90 °F to 160 °F).

- Silent Operation
- No Lint Traps required



Application

The Getinge GET250, GET350 and GET750 are designed to elevate blanket temperatures to a level which increases patient comfort.

Construction

- Stainless Steel exterior and white epoxy-coated interior
- Fully gasketed windowed doors with magnetic closure
- Four non-skid rubber feet are standard
- Right-hinged standard

Dimensions

| Model | Exterior (H x W x D) | Interior (H x W x D) | Net Weight |
|--------|---|--|-----------------|
| GET250 | 22"(559mm) x 18.5"(470mm) x 23.5" (597mm) | 15.25"(387mm) x 15.5"(394mm) x 18.75"(476mm) | 67 lbs (30 kg) |
| GET350 | 28"(711mm) x 18.5"(470mm) x 23.5"(597mm) | 21.25"(540mm) x 15.5"(394mm) x 18.75"(476mm) | 91 lbs (41 kg) |
| GET750 | 35.5"(902mm) x 24"(610mm) x 25.25"(641mm) | 28.25"(718mm) x 21"(533mm) x 21.75"(552mm) | 134 lbs (61 kg) |

Capacity

| Model | Bottles | Bags | Blankets |
|--------|---------|------|----------|
| GET250 | N/A | N/A | 4-6 |
| GET350 | N/A | N/A | 10-12 |
| GET750 | N/A | N/A | 20-25 |

Installation

Stationary or mobile warming cabinets

Number of doors

- Single door with window

Installation alternatives

- Bumper and Caster Assembly
- Caster and Plate Assembly
- Solid Door
- Left-hinged
- Pass-Through Design on GET750

Other installation options

- Manual lock
- Seismic Mounting brackets (specific installation requirements apply)
- Stacking Kit

Clearance Requirements

- 3" (76mm) from rear
- 1" (25mm) from top and sides
- 3/4" (19mm) from bottom

Safety features

- The controller stores all settings and resumes operation using these settings in the event of a power failure.
- A warming shut-off system, separate from the controller, prevents overheating.

Other features

- Access point at rear panel allows the addition of data logging or temperature management hardware.
- WarmSafe™ heating system incorporates a multiple zone warming technology (Patent No: US 8,217,316; US 8,581,152) that heats where and when it is needed.
- Operate in Celsius or Fahrenheit
- LED display
- Built-in speaker for audible feedback
- Actual temperature button
- GET350 and GET750 include a heated center shelf

Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries • Copyright 2019

Getinge AB or its subsidiaries or affiliates • All rights reserved • MCV00076065 REV B • EXP 02/21

Sales Office: Getinge USA Sales, LLC, 1777 East Henrietta Road, Rochester, NY 14623

Manufactured By: Entherrmics Medical Systems, W164 N9221 Water Street, Menomonee Falls, WI 53051

www.getinge.com

GETINGE 