

# Maquet Lyra

## A stable, flexible mobile OR table

Sometimes, a reliable OR table is all your hospital needs. The Maquet Lyra Mobile OR Table offers the stability and functionality needed across surgical disciplines, all for a comfortable price.



### The advantages of Maquet Lyra include:

#### **Excellent balance of price and performance**

Why pay more for features that you don't need? This cost-effective mobile OR table is flexible and easy to use, comfortably positioning patients for most types of surgical procedures.

#### **Exceptional stability and safety**

The robust base ensures stability and safety for patients weighing up to 360 kg (794 lbs). The Maquet Lyra Table is capable of a broad range of angles and tilts for extraordinary cross-disciplinary flexibility.

#### **Versatility across surgical disciplines**

By combining universal usage features with a broad range of accessories, the Maquet Lyra Mobile OR Table is a practical table for all-around use. Maquet Lyra features the same interfaces as past Getinge Mobile OR tables and the compatibility with existing Getinge accessories ensures seamless integration into your surgical suite.

#### **Reliability and user-friendly operation:**

The special design of the castors allows the table to be used immediately, without being lowered onto a baseplate. The free-wheeling castors give the flexibility to move in both straight lines or laterally. Fewer electrical components simplify maintenance for cost-effective long-term use.

#### **Longitudinal shift improves access for better outcomes**

The manual longitudinal shift improves intraoperative diagnostics and treatment. It delivers both increased access for the surgeon, and simultaneously improves accessibility for the C-arm. 360° imaging is possible when paired with carbon fiber accessories.

## Construction features

- Stable base construction with four double swivel castors for easy movement and maneuvering (base can be locked via foot pedal).
- Base cover made of Sheet Molding Composite (SMC).
- The override panel is connected to a separate circuit in order to ensure that the table can be safely controlled even in the event of a major software malfunction.
- Identical interfaces on normal and on reverse side.

## Electrical specifications

- Specially designed rechargeable batteries, with a single charge lasting about one week in the operating room.
- Electronic charge monitoring with visual indicators.
- Batteries recharged from the mains supply.
- AC100/110 - 120/127/200/220/230 - 240 V. (adjustable), 50/60 Hz, over mains cable.
- Safety class II, type B; the enclosure leakage current meets the requirements of the patient leakage current for type B conditions as per EN 60601-1.

### Characteristics Maquet Lyra 1118.10XX

Net weight  
(without accessories) 190 kg / 419 lbs

Length of OR table  
(without accessories) 745 mm / 29"

Maximum overall  
load (patient plus  
accessories) 360 kg / 794 lbs with restrictions  
180 kg / 397 lbs without  
restrictions

Trendelenburg/  
Reverse Trendelenburg 25° / 25°

Lateral tilt 15°

Lower back plate +75° / -40°

Leg plate (manual) -90° (down)

Manual longitudinal shift 310 mm / 12"



Getinge is a registered trademark of Getinge AB, its subsidiaries, or affiliates in the United States or other countries. ·  
Copyright 2020 Getinge AB or its subsidiaries or affiliates · All rights reserved.

**Manufactured by: Maquet (Suzhou) Co., Ltd.** · No.158 Fangzhou Road, Suzhou, P.R. China 215024 · +86 (0) 512 62839880

**Getinge** · 45 Barbour Pond Drive · Wayne, NJ 07470

**Sales Office, Canada** · 90 Matheson Blvd West, Suite 300 · Mississauga, Ontario L5R 3R3 · Canada

[www.getinge.com](http://www.getinge.com)

**GETINGE** 

MCV00103609 REVA EXP 05/2022