




Mega IAB Catheter

7.5 Fr. and 8 Fr. IAB Technical Specifications

Technical Specifications

Mega 7.5 Fr. and 8 Fr. IAB Catheter with Insertion Kit

Part Number	Balloon				
	Nominal Volume (cc)	Nominal Length (mm)	Balloon Dia. at 2 psi (mm)	Insertable Length (mm)	Guidewire Size mm (in.)
0684-00-0294-01	30	178	16	723	0.64 mm (0.025")
0684-00-0295-01	40	229	16	723	0.64 mm (0.025")
0684-00-0296-01	50	258	17.4	723	0.64 mm (0.025")

Ordering Information

Product	Part Number
Mega 7.5 Fr. 30 cc Balloon Catheter, Insertion Kit and two StatLock® IABs	0684-00-0294-01
Mega 7.5 Fr. 30 cc Balloon Catheter, Insertion Kit and two StatLock® IABs, with APA*	0684-00-0294-02
Mega 7.5 Fr. 40 cc Balloon Catheter, Insertion Kit and two StatLock® IABs	0684-00-0295-01
Mega 7.5 Fr. 40 cc Balloon Catheter, Insertion Kit and two StatLock® IABs, with APA*	0684-00-0295-02
Mega 8 Fr. 50 cc Balloon Catheter, Insertion Kit and two StatLock® IABs	0684-00-0296-01
Mega 8 Fr. 50 cc Balloon Catheter, Insertion Kit and two StatLock® IABs, with APA*	0684-00-0296-02

*Arrow Pump Adapter

Accessories	Part Number
Mega 7.5 Fr. 30 cc/40 cc Insertion Kit only	0884-00-0019-21
Mega 8 Fr. 50 cc Insertion Kit only	0884-00-0019-17

Mega Insertion Kits include:

- One 15 cm (6") Reinforced Introducer Sheath with hemostasis valve
- One Introducer Dilator
- One Vessel Dilator
- One 1.3 mm (18 G, 0.050") Angiographic Needle
- One 0.64 mm (0.025") X 145 cm, 3 mm J, PTFE coated, Stainless Steel Guidewire
- One Luer Cap
- One Three-Way Stopcock
- Two 122 cm (4') Pressure Tubings
- One 168 cm (5'6") Catheter Extender Tubing (30 cc/40 cc ONLY)
- One 170 cm (5'7") Catheter Extender Tubing (50 cc ONLY)

168 cm (5'6") Catheter Extender Tubing for 30 cc/40 cc IABs (10 Pack)	0684-00-0275
170 cm (5'7") Catheter Extender Tubing for 50 cc IABs (10 Pack)	0684-00-0526
StatLock® for IAB (20 pack) – Catheter Stabilization Device	0684-00-0472
Arrow Pump Adapter for 25 cc – 40 cc IABs (5 Pack)	0684-00-0510-01
Arrow Pump Adapter for 50 cc IABs (5 Pack)	0684-00-0510-02

Configuration

Packaging

One sterile pouch, corner peel

Cycle Life


One procedure

IAB Shelf Life

3 years from date of manufacturing

For more product information, contact your local representative.

The right size for every height



MCA00003496 REV B · Getinge, and **GETINGE**  are trademarks or registered trademarks of Getinge AB, its subsidiaries, or affiliates in the United States or other countries · StatLock is a trademark or registered trademark of C.R. Bard, Inc. or an affiliate · Copyright 2025 Datascope Corp. · All rights reserved ·  Refer to Instructions for Use for current indications, warnings, contraindications and precautions · 04/2025

 Getinge · 45 Barbour Pond Drive · Wayne, NJ 07470 · USA
 Datascope Corp. · 15 Law Drive · Fairfield, NJ 07004 · USA

www.getinge.com

GETINGE 