



# Getinge Smart Trolley range

Ergonomic flow of sterile goods



A trolley with Smart Options has integrated controls in the handle to enable safe use and maneuvering without letting go.

# Smart moves

## Different people, different positions

The Getinge Smart concept – a whole new level of convenience

### **An ergonomic concept**

Loading and distribution systems are an important part of the workflow in operating rooms and sterile supply units. Getinge Smart is a line of carefully designed trolleys that takes loading, unloading and distribution of instruments to a whole new level, both in terms of ergonomics and function.

The trolley range is based on an ergonomic concept for efficient loading and distribution of goods with the highest hygiene and usability. Getinge Smart means top-of-the-line functionality in

every detail, from the choice of materials and surfaces to maneuverability, cleanability and ergonomics.

Getinge Smart Trolleys can be equipped with a semi-automatic loading function to reduce manual handling of goods when loading or unloading. The trolleys are a joy to handle. The Power Drive function is an option enabling the user to maneuver the trolley to the desired position with little effort. These optional features assure good ergonomics and give your staff the best possible working conditions.

# Smart Loading and Unloading Trolleys

– ergonomic in every detail

Ergonomics is about people being able to work without exertion, strain or injury. It's about moving heavy loads easily. It's about being flexible and understanding that people differ in height and strength – and in the tasks they perform. Our loading and unloading trolleys\* are based on a Smart ergonomic concept that covers every detail of loading and distribution.

Smart Loading Trolleys without Smart Options can be washed and disinfected in Getinge 9100-series Cart Washer-Disinfectors (except trolley for Getinge HS33 Sterilizers).

## Smart fixed-height trolleys

- Easy to maneuver
- Easy to clean
- Ergonomic



Smart fixed-height trolley for washer-disinfectors



Smart fixed-height trolley for sterilizers

\* Containers, shelf racks and basket feeders need to be ordered separately.

# Smart Options

A trolley with Smart Options has integrated controls in the handle to enable safe use and maneuvering without letting go.

## **Semi-automatic loading (SAL)**

For automatically loading and unloading sterilizer racks, activated by the control integrated in the handle.

## **Smart adjustable-height trolleys (AHT)**

Adjust to ergonomic working heights to reduce strain for all users.

## **Power drive (PD)**

Significantly reduced strain by enabling moving and maneuvering of trolleys to the desired position with little effort.



Smart adjustable-height trolley for washer-disinfectors



Smart adjustable-height trolley for sterilizers

# Smart Loading and Unloading Trolleys

## Technical information

### Loading trolley for rack – Getinge HS33 Sterilizers

| Sterilizer | Item no.  | Product description | Max. load      |
|------------|-----------|---------------------|----------------|
| HS33-60L   | 564710370 | Smart FHT HS33-60L  | 30 kg / 66 lbs |
| HS33-100L  | 564710372 | Smart FHT HS33-100L | 30 kg / 66 lbs |



### Loading trolley for rack – Getinge 522/533-HC Sterilizers

| Sterilizer | Item no.  | Product description | Max. load        |
|------------|-----------|---------------------|------------------|
| 522/533-HC | 564710471 | LT FIX FOR 500      | 130 kg / 287 lbs |



### Loading trolley for shelf rack – Getinge 66, 67, 633HC and Solsus 66 (ST001) Sterilizers

#### Full-length load carrier

| Sterilizer                               | Item no.   | Product description            | Max. load        | Length load carrier |
|--|------------|--------------------------------|------------------|---------------------|
| HS/GE 6606, 633-26"                      | 6020720670 | Smart FHT 66/67x06/13HL/633-07 | 200 kg / 440 lbs | 640 mm / 25.2 in    |
| HS/GE 6610, 633-39"<br>GSS6710, ST001-10 | 6020720672 | Smart FHT 66/67x10/20HL/633-10 | 200 kg / 440 lbs | 990 mm / 38.97 in   |
| HS/GE 6613, 633-51"<br>GSS6713, ST001-13 | 6020720674 | Smart FHT 66/67x13/633-13      | 200 kg / 440 lbs | 1,290 mm / 50.78 in |
| HS/GE 6617,<br>GSS6717, ST001-17         | 6020720676 | Smart FHT 66/67x17             | 200 kg / 440 lbs | 1,680 mm / 6.14 in  |



#### Half-length load carrier – Two load carriers for full sterilizer load

| Sterilizer                               | Item no.   | Product description            | Max. load        | Length load carrier |
|--|------------|--------------------------------|------------------|---------------------|
| HS/GE 6613, 633-51"<br>GSS6713, ST001-13 | 6020720670 | Smart FHT 66/67x06/13HL/633-07 | 200 kg / 440 lbs | 640 mm / 25.2 in    |
| HS/GE 6617,<br>GSS6717, ST001-17         | 6020720672 | Smart FHT 66/67x10/20HL/633-10 | 200 kg / 440 lbs | 840 mm / 33.07 in   |
| HS/GE 6620,<br>GSS6720, ST001-20         | 6020720672 | Smart FHT 66/67x10/20HL/633-10 | 200 kg / 440 lbs | 840 mm / 33.07 in   |

Getinge Smart Loading and Unloading Trolleys\* give you the possibility to choose a solution that really meets your needs in every situation.



## Loading trolley for shelf rack – Getinge 66, 67, 633HC and Solsus 66 (ST001) Sterilizers

| With Smart Option full-length load carrier |            |                                |                  |                     |     | Smart Options |    |  |
|--|------------|--------------------------------|------------------|---------------------|-----|---------------|----|--|
| Sterilizer                                 | Item no.   | Product description            | Max. load        | Length load carrier | AHT | SAL           | PD |  |
| HS/GE 6610 GSS6710 ST001-10                | 6020721472 | Smart FHT SAL 66/67x10, 633-10 | 200 kg / 440 lbs | 990 mm / 38.97 in   | -   | •             | -  |  |
| HS/GE 6613 GSS6713 ST001-13                | 6020721474 | Smart FHT SAL 66/67x13, 633-13 | 200 kg / 440 lbs | 1,290 mm / 50.78 in | -   | •             | -  |  |
| HS/GE 6610 GSS6710 ST001-10                | 6020723072 | Smart AHT 66/67/69x10, 633-10  | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | -             | -  |  |
| HS/GE 6613 GSS6713 ST001-13                | 6020723074 | Smart AHT 66/67/69x13, 633-13  | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •   | -             | -  |  |
| HS/GE 6617 GSS6717 ST001-17                | 6020723076 | Smart AHT 66/67x17             | 200 kg / 440 lbs | 1,680 mm / 66.14 in | •   | -             | -  |  |
| HS/GE 6610 GSS6710 ST001-10                | 6020724072 | Smart AHT SAL 6x10/633-10      | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | •             | -  |  |
| HS/GE 6613 GSS6713 ST001-13                | 6020724074 | Smart AHT SAL 6x13/633-13      | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •   | •             | -  |  |
| HS/GE 6610 GSS6710 ST001-10                | 6020724572 | Smart AHT PD 6x10/633-10       | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | -             | •  |  |
| HS/GE 6613 GSS6713 ST001-13                | 6020724574 | Smart AHT PD 6x13/633-13       | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •   | -             | •  |  |
| HS/GE 6610 GSS6710 ST001-10                | 6020725072 | Smart AHT SAL PD 6x10/633-10   | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | •             | •  |  |
| HS/GE 6613 GSS6713 ST001-13                | 6020725074 | Smart AHT SAL PD 6x13/633-13   | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •   | •             | •  |  |

| With Smart Option half-length load carrier – Two load carriers for full sterilizer load |            |                               |                  |                     |     | Smart Options |    |  |
|---|------------|-------------------------------|------------------|---------------------|-----|---------------|----|--|
| Sterilizer  | Item no.   | Product description           | Max. load        | Length load carrier | AHT | SAL           | PD |  |
| HS/GE 6617, GSS6717 ST001-17  | 6020723072 | Smart AHT 66/67/69x10, 633-10 | 200 kg / 440 lbs | 840 mm / 33.07 in   | •   | -             | -  |  |
| HS/GE 6617, GSS6717 ST001-17  | 6020724572 | Smart AHT PD 6x10/633-10      | 200 kg / 440 lbs | 840 mm / 33.07 in   | •   | -             | •  |  |
| HS/GE 6620, GSS6720 ST001-20  | 6020723072 | Smart AHT 66/67/69x10, 633-10 | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | -             | -  |  |
| HS/GE 6620, GSS6720 ST001-20  | 6020724572 | Smart AHT PD 6x10/633-10      | 200 kg / 440 lbs | 990 mm / 38.97 in   | •   | -             | •  |  |

| Only for docking to automatic loading system, for shorter footprint |            |                             |                  |                     |     | Smart Options |    |  |
|---|------------|-----------------------------|------------------|---------------------|-----|---------------|----|--|
| Sterilizer  | Item no.   | Product description         | Max. load        | Length load carrier | AHT | SAL           | PD |  |
| HS/GE 6617 GSS6717  | 6020723176 | Smart AHT 66/67x17 HL UN-/L | 200 kg / 440 lbs | 840 mm / 33.07 in   | •   | -             | -  |  |

\* Containers, shelf racks and basket feeders need to be ordered separately.

## Loading trolley for basket feeder – Getinge 66, 67 and Solsus 66 (ST001) Sterilizers

The loading trolley together with a basket feeder enable manual loading of baskets into a sterilizer

| Sterilizer                      | Item no.   | Product description            | Max. load         | Basket feeder |
|---------------------------------|------------|--------------------------------|-------------------|---------------|
| HS/GE 6610, GSS6710<br>ST001-10 | 6020720672 | Smart FHT 66/67x10/20HL/633-10 | 125 kg / 275 lbs* | 564613270CE   |
| HS/GE 6613, GSS6713<br>ST001-13 | 6020720674 | Smart FHT 66/67x13/633-13      | 125 kg / 275 lbs* | 564613470CE   |
| HS/GE 6617, GSS6717<br>ST001-17 | 6020720676 | Smart FHT 66/67x17             | 125 kg / 275 lbs* | 564613670CE   |



### With Smart Option

| Sterilizer                      | Item no.   | Product description           | Max. load         | Basket feeder | Smart Options |     |    |
|---------------------------------|------------|-------------------------------|-------------------|---------------|---------------|-----|----|
|                                 |            |                               |                   |               | AHT           | SAL | PD |
| HS/GE 6610, GSS6710<br>ST001-10 | 6020723072 | Smart AHT 66/67/69x10, 633-10 | 125 kg / 275 lbs* | 564613270CE   | ●             | -   | -  |
| HS/GE 6610, GSS6710<br>ST001-10 | 6020724572 | Smart AHT PD 6x10/633-10      | 125 kg / 275 lbs* | 564613270CE   | ●             | -   | ●  |
| HS/GE 6613, GSS6713<br>ST001-13 | 6020723074 | Smart AHT 66/67/69x13, 633-13 | 125 kg / 275 lbs* | 564613470CE   | ●             | -   | -  |
| HS/GE 6613, GSS6713<br>ST001-13 | 6020724574 | Smart AHT PD 6x13/633-13      | 125 kg / 275 lbs* | 564613470CE   | ●             | -   | ●  |
| HS/GE 6617, GSS6717<br>ST001-17 | 6020723076 | Smart AHT 66/67x17            | 125 kg / 275 lbs* | 564613670CE   | ●             | -   | -  |

\*The capacity of the basket is limited to 125 kg / 275 lbs.

## Loading trolley for automatic basket loader / unloader – Getinge 66 and 67 Sterilizers

The loading trolley is used to load / unload an automatic loader and unloader for baskets.  
Not possible to dock directly to sterilizer.



### With Smart Option

| Sterilizer          | Item no.   | Product description             | Max. load        | Smart Options |     |    |
|---------------------|------------|---------------------------------|------------------|---------------|-----|----|
|                     |            |                                 |                  | AHT           | SAL | PD |
| HS/GE 6610, GSS6710 | 6020727072 | Smart AHT basket 66/67x10 UN-/L | 200 kg / 440 lbs | ●             | -   | -  |
| HS/GE 6613, GSS6713 | 6020727074 | Smart AHT basket 66/67x13 UN-/L | 200 kg / 440 lbs | ●             | -   | -  |



## Loading trolley for shelf racks – Getinge 69 Sterilizers

### Full-length load carrier

| Sterilizer | Item no.   | Product description | Max. load        | Length load carrier |
|------------|------------|---------------------|------------------|---------------------|
| HS/GE 6910 | 6020720872 | Smart FHT 6910      | 200 kg / 440 lbs | 990 mm / 38.97 in   |
| HS/GE 6913 | 6020720874 | Smart FHT 6913      | 200 kg / 440 lbs | 1,290 mm / 50.78 in |
| HS/GE 6915 | 6020720875 | Smart FHT 6915      | 210 kg / 462 lbs | 1,530 mm / 60.23 in |



## Loading trolley for shelf racks – Getinge 69 Sterilizers

### With Smart Option full-length load carrier

| Sterilizer | Item no.   | Product description           | Max. load        | Length load carrier | Smart Options |     |    |
|------------|------------|-------------------------------|------------------|---------------------|---------------|-----|----|
|            |            |                               |                  |                     | AHT           | SAL | PD |
| HS/GE 6910 | 6020721672 | Smart FHT SAL 6910            | 200 kg / 440 lbs | 990 mm / 38.97 in   | -             | •   | -  |
| HS/GE 6913 | 6020721674 | Smart FHT SAL 6913            | 200 kg / 440 lbs | 1,290 mm / 50.78 in | -             | •   | -  |
| HS/GE 6910 | 6020723072 | Smart AHT 66/67/69x10, 633-10 | 200 kg / 440 lbs | 990 mm / 38.97 in   | •             | -   | -  |
| HS/GE 6913 | 6020723074 | Smart AHT 66/67/69x13, 633-13 | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •             | -   | -  |
| GE 6915    | 6020723075 | Smart AHT 6915                | 210 kg / 462 lbs | 1,530 mm / 60.23 in | •             | -   | -  |
| HS/GE 6910 | 6020724072 | Smart AHT SAL 6x10/633-10     | 200 kg / 440 lbs | 990 mm / 38.97 in   | •             | •   | -  |
| HS/GE 6913 | 6020724074 | Smart AHT SAL 6x13/633-13     | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •             | •   | -  |
| HS/GE 6610 | 6020724572 | Smart AHT PD 6x10/633-10      | 200 kg / 440 lbs | 990 mm / 38.97 in   | •             | -   | •  |
| HS/GE 6613 | 6020724574 | Smart AHT PD 6x13/633-13      | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •             | -   | •  |
| GE 6915    | 6020724575 | Smart AHT PD 6915             | 210 kg / 462 lbs | 1,530 mm / 60.23 in | •             | -   | •  |
| HS/GE 6910 | 6020725072 | Smart AHT SAL PD 6x10/633-10  | 200 kg / 440 lbs | 990 mm / 38.97 in   | •             | •   | •  |
| HS/GE 6913 | 6020725074 | Smart AHT SAL PD 6x13/633-13  | 200 kg / 440 lbs | 1,290 mm / 50.78 in | •             | •   | •  |
| GE 6915    | 6020725075 | Smart AHT SAL PD 6915         | 210 kg / 462 lbs | 1,530 mm / 60.23 in | •             | •   | •  |



### With Smart Option half-length load carrier – Two load carriers for full sterilizer load

| Sterilizer | Item no.   | Product description           | Max. load        | Length load carrier | Options |     |    |
|------------|------------|-------------------------------|------------------|---------------------|---------|-----|----|
|            |            |                               |                  |                     | AHT     | SAL | PD |
| HS/GE 6917 | 6020723072 | Smart AHT 66/67/69x10, 633-10 | 200 kg / 440 lbs | 840 mm / 33.07 in   | •       | -   | -  |

## Loading trolley for shelf racks – Getinge 69 GEB Sterilizers

Load carrier Getinge Shelf Rack, loading height 1,100 mm / 40.30 in.  
Possible to dock to sterilizer.



### With Smart Option full-length load carrier

### Smart Options

| Sterilizer | Item no.   | Product description | Max. load        | Load length carrier | AHT | SAL | PD |
|------------|------------|---------------------|------------------|---------------------|-----|-----|----|
| GEB 6910   | 6020726072 | Smart AHT GEB 6910  | 200 kg / 440 lbs | 990 mm / 38.97 in   | ●   | -   | -  |
| GEB 6913   | 6020726074 | Smart AHT GEB 6913  | 200 kg / 440 lbs | 1,290 mm / 50.78 in | ●   | -   | -  |
| GEB 6915   | 6020726075 | Smart AHT GEB 6915  | 210 kg / 462 lbs | 990 mm / 38.97 in   | ●   | -   | -  |

## Loading trolley for shelf racks – Getinge 610 Sterilizers

Load carrier Getinge Shelf Rack. Possible to dock to sterilizer.

### Full-length load carrier

| Sterilizer | Item no.   | Product description             | Max. load        | Length load carrier |
|------------|------------|---------------------------------|------------------|---------------------|
| GSS61010   | 6023373272 | Smart Load Trolley FHT 61010/21 | 265 kg / 585 lbs | 1,000 mm / 39.37 in |
| GSS61014   | 6023373274 | Smart Load Trolley FHT 61014    | 265 kg / 585 lbs | 1,360 mm / 53.54 in |
| GSS61015   | 6023373275 | Smart Load Trolley FHT 61015    | 265 kg / 585 lbs | 1,500 mm / 59.05 in |




### Half-length load carrier – Two load carriers for full sterilizer load

| Sterilizer | Item no.   | Product description               | Max. load        | Length load carrier |
|------------|------------|-----------------------------------|------------------|---------------------|
| GSS61014   | 6023373276 | Smart Load Trolley FHT 61014/15HL | 265 kg / 585 lbs | 640 mm / 25.2 in    |
| GSS61015   | 6023373276 | Smart Load Trolley FHT 61014/15HL | 265 kg / 585 lbs | 710 mm / 27.9 in    |
| GSS61021   | 6023373272 | Smart Load Trolley FHT 61010/21   | 265 kg / 585 lbs | 1,000 mm / 39.37 in |

## Loading trolleys for washer-disinfectors

Load carrier Getinge Wash Cart.\* Possible to dock to 8666/8668/86T/88T/CM320 Washer-Disinfectors and loading equipment for Automation System 2.0.



Smart Options

### With Smart Option single-end loading

| Washer-Disinfector      | Item no.   | Product description                     | Max. load      | AHT | SAL | PD |
|-------------------------|------------|---|----------------|-----|-----|----|
| 8666/8668/86T/88T/CM320 | 6020720271 | Smart AHT WD<br>8666/8668/86T/88T/CM320 | 165 kg/363 lbs | •   | -   | -  |

### Single-end loading

| Washer-Disinfector      | Item no.   | Product description                     | Max. load      |
|-------------------------|------------|---|----------------|
| 8666/8668/86T/88T/CM320 | 6020720171 | Smart FHT WD<br>8666/8668/86T/88T/CM320 | 165 kg/363 lbs |

### Double-end loading

| Washer-Disinfector      | Item no.   | Product description                        | Max. load      |
|-------------------------|------------|--|----------------|
| 8666/8668/86T/88T/CM320 | 6020720172 | Smart FHT WD DH<br>8666/8668/86T/88T/CM320 | 165 kg/363 lbs |



\*Wash cart needs to be ordered separately.

# Smart Distribution Trolleys

– for ergonomic transport efficiency

With Getinge Smart Distribution Trolleys you always get safe, efficient distribution. The hard-cover model protects sterile goods from dust and damage. The open trolley model gives you a mobile unit for modular baskets and containers.

## Smart Distribution Trolley Hard-Cover (DTHC)

In this trolley, instrument containers and baskets can be grouped together with utensils, linen, disposables and anesthesia material for transport to an internal OR department or to an external hospital or storage facility. The hard cover protects the sterile goods from dust and damage during transport.

- Simple and safe transport
- Ergonomic
- Can be washed and disinfected in Getinge Cart Washer-Disinfector



Getinge Smart Distribution Trolleys are designed for the effective, ergonomic transport of baskets and/or containers. They are used in packing areas in sterile zones, and for transport to and from soiled zones throughout the hospital.



The doors can be opened 260° and fixed in the open position – very convenient when loading and unloading. The door handle can be fitted with an optional lock. Shelf supports can be moved, which maximizes utilization of the trolley.



A maximum of five shelves can be used. Three are supplied as standard. Additional shelves are available (optional) for extra storage space. Each shelf is supplied with two dividers to keep goods in place during transport.



To facilitate maneuvering, the trolley is equipped with two lockable swivel casters and two rigid casters. The caster diameter is 160 mm. An optional pin and tow bar is available to enable towing.

### Smart Distribution Trolley Hard-Cover (DTHC)

| Item no.  | Capacity  | Options (Item no.)   | Max. load      | Trolley dimensions (L x W x H)            |
|-----------|---|--|----------------|---|
| 564719070 | 6 or 8 SPRI/ISO high*<br>or 15 DIN StU ½ size** | 564719071 wire shelf incl. 2 dividers<br>564719075 pin and tow bar | 145 kg/319 lbs | 1,286x699x1,502 mm / 50.63x27.52x59.13 in |

\* SPRI high: 585x395x195 mm / 23.03x15.55x7x67 in. ISO high: 600x400x195 mm / 23.62x15.75x7x67 in (for 8 SPRI/ISO, order one extra shelf)

\*\* DIN StU ½ size: 575x280x135 mm / 22.64x11.02x5.31 in

# Smart Distribution Trolleys Open



## Smart Distribution Trolley Open (DTO)

This trolley is a mobile storage unit for modular baskets and containers. The open design makes viewing and picking of individual packs very easy. Shelves and baskets are easily secured for safe transport.

- Easy to load and unload
- Saves storage space
- Silent and easy to maneuver
- Can be washed and disinfected in Getinge Cart Washer-Disinfector



The flexible placement of support guides and shelves maximizes trolley utilization.



Trolleys can be loaded and unloaded from both sides. Shelves and baskets can be pulled out halfway for easier access.



The high-quality swivel casters (two of which are lockable) enable easy moving and maneuvering even when the trolley is fully loaded.




**DTO double (shelves)**  
Item no.: 564718970



**DTO single (mixed)**  
Item no.: 564718370



**DTO single SPRI (5 high)**  
Item no.: 564718170  
(Baskets not included)



**DTO single (shelves)**  
Item no.: 564718570  
(Container not included)

**Getinge Smart Distribution Trolleys Open are available in single or double versions of each model:**

- Three models configured for standard loads. Delivered with guides for baskets, shelves for containers or mixed
- One base model that can be individually configured. For the base model, guides are specified and ordered separately for customers' individual configuration of a mixed load of shelves and /or baskets.

**Technical data – Smart Distribution Trolley Open (DTO)**

**DTO single**

| Type of load            | Item no.  | Product description                 | Max. load       | Dimensions (LxWxH)                     |
|-------------------------|-----------|-------------------------------------|-----------------|--|
| Baskets                 | 564718170 | DTO single SPRI, 10 low or 5 high   | 85 kg / 187 lbs | 736x429x1577 mm / 28.98x16.89x62.09 in |
| Mixed                   | 564718370 | DTO single, 2 SPRI high / 3 shelves | 85 kg / 187 lbs | 736x429x1577 mm / 28.98x16.89x62.09 in |
| Shelves                 | 564718570 | DTO single, 5 shelves               | 85 kg / 187 lbs | 736x429x1577 mm / 28.98x16.89x62.09 in |
| Individually configured | 564718070 | DTO single base without guides      | 85 kg / 187 lbs | 736x429x1577 mm / 28.98x16.89x62.09 in |

**DTO double**

| Type of load            | Item no.  | Product description                 | Max. load        | Dimensions (LxWxH)                        |
|-------------------------|-----------|-------------------------------------|------------------|---|
| Baskets                 | 564718770 | DTO double, SPRI, 20 low or 10 high | 170 kg / 374 lbs | 1,416x429x1,577 mm / 55.75x16.89x62.09 in |
| Mixed                   | 564718870 | DTO double, 4 SPRI high / 6 shelves | 170 kg / 374 lbs | 1,416x429x1,577 mm / 55.75x16.89x62.09 in |
| Shelves                 | 564718970 | DTO double, 10 shelves              | 170 kg / 374 lbs | 1,416x429x1,577 mm / 55.75x16.89x62.09 in |
| Individually configured | 564718670 | DTO double base without guides      | 170 kg / 374 lbs | 1,416x429x1,577 mm / 55.75x16.89x62.09 in |

**Options for single and double**

| Item no.   | Product description                          | Max. load      | Dimensions (LxWxH)                |
|------------|--|----------------|-----------------------------------|
| 564718171  | Pair of guides SPRI                          | 18 kg / 39 lbs | 30x425x50 mm / 1.18x16.73x1.99 in |
| 564718771  | Pair of guides ISO                           | 18 kg / 39 lbs | 30x425x50 mm / 1.18x16.73x1.99 in |
| 564718571  | Pair of guides and shelf                     | 18 kg / 39 lbs | 30x425x50 mm / 1.18x16.73x1.99 in |
| 564718470  | External transport lock                      |                |                                   |
| 560181301  | Rack for disposable covers                   |                |                                   |
| 560125070  | Disposable transparent covers, 50 pcs / roll |                |                                   |
| 6003000052 | Reusable cover for DTO single                |                |                                   |
| 6003000053 | Reusable cover for DTO double                |                |                                   |



With a firm belief that every person and community should have access to the best possible care, Getinge provides hospitals and life science institutions with products and solutions aiming to improve clinical results and optimize workflows. The offering includes products and solutions for intensive care, cardiovascular procedures, operating rooms, sterile reprocessing and life science. Getinge employs over 10,000 people worldwide and the products are sold in more than 135 countries.

This document is intended to provide information to an international audience outside of the US.

**Getinge Sterilization AB** · PO Box 69 · SE-305 05 Getinge · Sweden · +46 10 335 00 00 · [info@getinge.com](mailto:info@getinge.com)

[www.getinge.com](http://www.getinge.com)