

# Wash-Checks Pro Protein Detection Kit

## WC 112

Wash-Checks Pro is a quick and easy test used to detect residual protein on surfaces of surgical instruments and washing equipment.



### Key Features

- 3 inch swab (76 mm)
- The higher level of contamination of protein the darker and faster the blue result
- Accuracy of 1 µg in less than 10 seconds
- No incubation required

### Performance Characteristics

This protein test can detect protein residues within 1 µg sensitivity and 10 seconds. When protein is present, the solution will turn from brown to blue. This test provides a semi-quantitative measure of surface cleanliness; the higher the level of contamination of protein (bioburden), the darker and faster the color changes to blue.

### Storage

Storage conditions are printed on the product packaging. Keep away from light and heat.

### Shelf Life

2 years from manufacture of reagent when refrigerated and 6 months when not refrigerated.

### Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors

#### SteriTec Part No.

WC 112, Wash-Checks Pro

#### Packaging

2 kits of 25 tests & 25 swabs