

# Wash-Checks Pro Protein Detection Kit 2.5m Swab WC 118

Wash-Checks Pro is a quick and easy test used to detect residual protein in hollow lumen devices and endoscopes.



## Key Features

- 2.5m swab
- The higher level of contamination of protein the darker and faster the blue result
- Accuracy of 1 µg in less than 10 seconds
- No incubation required

## Performance Characteristics

This protein test can detect protein residues within 1 µg sensitivity and 10 seconds. When protein is present, the solution will turn from brown to blue. This test provides a semi-quantitative measure of surface cleanliness; the higher the level of contamination of protein (bioburden), the darker and faster the color changes to blue.

## Storage

Storage conditions are printed on the product packaging. Keep away from light and heat.

## Shelf Life

2 years from manufacture of reagent when refrigerated and 6 months when not refrigerated.

## Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors

### SteriTec Part No.

WC 118, Wash-Checks Pro

### Packaging

2 boxes of 10 tests & 10 swabs

FTP-350-2 Rev. B



SteriTec Products Mfg. Co., Inc.,  
Member of the Getinge Group  
74 Inverness Drive East  
Englewood, CO 80112 USA

[www.getinge.com/steritec](http://www.getinge.com/steritec)