

# Trimano Fortis

## Improving ergonomics

For years, the Trimano 3D Support Arm has offered a helping hand in the operating theater. Getinge continues this success story with the Trimano Fortis.



Trimano Fortis brings strength, flexibility and ergonomic advantages to the operating theater. It is a further development of the Trimano 3D support arm. Trimano Fortis is more robust than its predecessor, with significant improvements to load capacity\* and handling. Trimano Fortis builds upon the popular features of the widely used Trimano 3D Support Arm. It also offers ergonomic benefits for clinicians. The force required to release the handle has been reduced, making repositioning easier.

\* max. proportional weight of patient arm of 250 kg

### Benefits at a glance:

- **Easy to use:** Trimano Fortis requires no power supply or additional connection
- **Work relief:** thanks to its versatile applications the Trimano improves working conditions the surgical team
- **Extremely movable:** the "Click and Move" system enables free movement in all directions while securing a stable fixation when the handle is released
- **Multidisciplinary:** due to a large range of accessories, the Trimano is usable in many surgical disciplines
- **Compatibility:** can be used with common OR tables with a side rail
- **Maximize your investment:** backwards compatible with predecessor Trimano (with exceptions), meaning existing accessories can be used

# An extra-hand in the operating room

Holding and lifting body extremities during orthopedic interventions is a huge effort for surgical assistants. This task is performed by Trimano Fortis easily and with more stability.

The Trimano Fortis can be connected directly to the OR table accessory rail and is ready for immediate use. No need for any other connections. It can be moved in any direction and securely fixed at any position.

Trimano's features – such as multidisciplinary use, flexible positioning options, increased load capacity, and comfortable handling – make it an essential tool in the operating theater. And as Trimano Fortis is mostly compatible with the previous version, you can continue using your existing equipment.

## Ordering numbers and compatible accessories

1002.30A0	Trimano Fortis
1002.31A0	Adapter Fortis
1002.32A0	Armboard Fortis, SFC
1002.43A0	Arm rest, lateral position
1002.43B0	Arm rest, beach chair
1002.52A0	Arm rest, disposable, sterile
1002.45A0	Countertraction post
1002.48A0	Finger fixation

Disposables are not listed, please refer to the **OR tables accessories catalog**.



Shoulder treatment in beach chair position



Wrist treatment in supine position



Plastic surgery in lateral position



Elbow treatment in prone position

This document is intended to provide information to an international audience outside of the US.  
Maquet GmbH · Kehler Str. 31 · 76437 Rastatt · Germany · +49 7222 932-0

[www.getinge.com](http://www.getinge.com)

GETINGE 