


XEN

Chemistries
for the ultimate
peace of mind

Enzymatic+



Scan for product
information
including Safety
Data Sheets

California Ingredients Disclosure

Product Codes: 6036002901, 6036003001, 6036003101

Ingredients	CAS Number	Purpose	Additional Information
Water	7732-18-5	Diluent	N/A
Triethanolamine	102-71-6	Builder	https://open.canada.ca/data/en/dataset/1d946396-cf9a-4fa1-8942-4541063bfba4
Organic acid	Withheld	Withheld	N/A
Surfactant	Withheld	Withheld	N/A
Polymer	Withheld	Withheld	N/A
Surfactant	Withheld	Withheld	N/A
Subtilisin	9014-01-1	Cleaning Agent	https://echa.europa.eu/information-on-chemicals/cl-inventory-database
Surfactant	Withheld	Withheld	N/A
Preservative	2634-33-5	Preservation	https://open.canada.ca/data/en/dataset/1d946396-cf9a-4fa1-8942-4541063bfba4
Corrosion inhibitor	Withheld	Withheld	N/A
Colorant	Withheld	Withheld	N/A
Amylase	9000-90-2	Cleaning Agent	https://echa.europa.eu/information-on-chemicals/cl-inventory-database
Lipase	9001-62-1	Cleaning Agent	https://echa.europa.eu/information-on-chemicals/cl-inventory-database
Cellulase	9012-54-8	Cleaning Agent	https://echa.europa.eu/information-on-chemicals/cl-inventory-database



Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries.
Copyright 2024 Getinge AB or its subsidiaries or affiliates. • All rights reserved.

Manufacturer: Quadralene Ltd, Unit 3, St. Modwen Park, Andressey Way, Derby, DE21 6YH, United Kingdom.

Sales Office, US: 1 Geoffrey Way, Wayne, NJ 07470, USA

Sales Office, Canada: 90 Matheson Blvd West, Suite 300, Mississauga, Ontario L5R 3R3, Canada

www.getinge.com

GETINGE 

MCV00111532 REVA